

BioPhase 8800 system: Advancing capillary electrophoresis with native fluorescence detection

SCIEX CE team, 2025

Dr. Fang Wang, Sr. Product Manager, SCIEX, US

BioPhase 8800 system with NF

Year : **2025**

TAKE **CHARGE**



*Maintain all existing capabilities on the BioPhase 8800 system.
Benefit from the addition of native fluorescence detection:*

- 1. Protein detection with improved sensitivity compared to traditional UV detection.*
- 2. Faster, more confident peak integration than traditional UV due to significantly reduced noise and flattened baseline.*
- 3. No need for protein-dye conjugation during sample preparation, speeding up sample preparation.*



Features - hardware

- Factory-built cartridge with push-in load: Reduce user error
- 8 capillaries parallel processing: Speed up your sample analysis
- Built-in electrodes, interface block, and aperture: Reduce maintenance
- Smart chip with SN, usage tracking: Improve traceability
- Precise temperature capillary control: Maximize data reproducibility
- Plate-based sample/reagent: Automate processes
- Software-driven detector selection: Simplify detector changes and enable a single sequence with multiple detector types
- Reduced stray light and higher transmission: Gain higher signal intensity



Detector Type

☐ UV	Wavelength	220	nm
☐ Wait			
☐ LIF	Emission Wavelength	520	nm
☐ Wait	PMT Gain	100	
☒ NFD	Emission Wavelength	350	nm
☒ Wait	Application	Other	
	PMT Gain	100	






Figure 3. Schematic diagram of the stray light.

Features - software


The user can choose to set up either in a local or a network configuration




Intuitive touchscreen interface



Drag and drop method programming



BioPhase 8800 Driver for Empower software & documentation + OQ3



Native fluorescence of proteins: A new way of detection

The native fluorescence of proteins originates from residues of three aromatic amino acids:
tryptophan , tyrosine, and phenylalanine

Tryptophan : $\lambda_{ex}/\lambda_{em}=280/350$ nm

-- High emission intensity, the primary native fluorescence source in proteins

Tyrosine : $\lambda_{ex}/\lambda_{em}=284/304$ nm & $\lambda_{ex}/\lambda_{em}=220/304$ nm

-- High emission intensity, but quenched by its interaction with the peptide chain via energy transfer

Phenylalanine : $\lambda_{ex}/\lambda_{em}=250/280$ nm

-- Low emission intensity


No Need to rely on low-wavelength amide bond or fluorescent tag labeling for protein detection

A full-page background image of a male scientist with a long white beard and safety glasses, wearing a white lab coat and blue gloves. He is using a pipette to transfer liquid into a clear bottle. The background shows a laboratory environment with various equipment and papers on the wall.

CE-SDS analysis for proteins:

BioPhase 8800 system with NF: non-reduced CE-SDS

Traditional UV 220 nm detection




Protein labeling is not needed

Requires high sample concentration

Non-reduced CE-SDS

Sample: 1 mg/mL [1 mg/mL] non-reduced NIST IgG


NEW! NF 350 nm detection



Protein labeling is not needed

Improved S/N ratio & flat baseline
compared to UV detection

LIF 600 nm detection



LIF dye conjugation to protein is required

Improved S/N ratio > 10X & flat baseline
compared to UV detection


BioPhase 8800 system with NF: non-reduced purity for mAb

- The measured relative abundance of total fragments was comparable among the three detection approaches
- The measured relative abundance of aggregates with NFD is equivalent to that with LIF
- The aggregate peaks were not detectable with UV due to interference by the baseline anomalies


Sample Lot	Detection	Relative Abundance, %							
		L	H	HL	HH	HHL	Total	Intact mAb	Aggregates
1	UV	0.29	0.13	0.15	0.18	0.78	1.53	98.48	ND*
	NFD	0.24	0.14	0.15	0.24	1.00	1.77	97.84	0.40
	LIF	0.18	0.15	0.15	0.22	0.92	1.63	97.97	0.40
2	UV	0.29	0.14	0.11	0.17	1.07	1.78	98.22	ND*
	NFD	0.22	0.17	0.14	0.21	1.26	2.00	97.62	0.38
	LIF	0.24	0.14	0.14	0.25	1.09	1.86	97.77	0.37
3	UV	0.28	0.19	0.13	0.18	0.96	1.73	98.27	ND*
	NFD	0.21	0.20	0.15	0.22	1.11	1.88	97.74	0.38
	LIF	0.24	0.14	0.16	0.26	1.06	1.85	97.80	0.35

**ND – not detected. With UV absorbance detection, the aggregate peaks were not detected due to interference by the baseline anomalies.*

BioPhase 8800 system with NF: Quantitation of 0.01% spiked in lysozyme



NF 350 nm




Sample: non-reduced USP IgG at 1 mg/mL.

- NFD 350 nm: S/N RMS at 0.1% spike in > 1000, linear extrapolation S/N RMS at 0.01% spike in > 10.
- UV 220 nm: S/N RMS at 0.1% spike in ~ 20


BioPhase 8800 system with NF: Quantification of protein fragments

Traditional UV detection: USP IgG



NEW!

NF detection USP IgG




C

	UV detection		NF detection	
	CPA %	Average S/N RMS	CPA %	Average S/N RMS
USP IgG-HC	0.53	78	0.44	185
NIST mAb- HL	0.25	20	0.19	76
NIST mAb-HH	0.22	19	0.21	83

BioPhase 8800 system with NF: reduced CE-SDS

Traditional UV 220 nm detection



Protein labeling is not needed


Requires high sample concentration

Reduced CE-SDS

Sample: 1000 µg/mL [1 mg/mL] reduced NIST IgG

NEW!


NFD 350 nm detection



Protein labeling is not needed
and

Improved sensitivity & flat baseline
compared to UV detection

LIF-600 nm detection



LIF dye conjugation to protein is required

Improved sensitivity >10X & flat baseline
compared to UV detection

BioPhase 8800 system with NF: glycan occupancy


Sample Lot	Glycosylation Site Occupancy on Heavy Chain, %		
	UV Absorbance	NFD	LIF
1	99.30	99.33	99.33
2	99.31	99.32	99.32
3	99.31	99.33	99.33

The determined glycosylation site occupancy was consistent, in the range of 99.30–99.33%, using the three detection approaches.

NIST mAb
Glycosylation site occupancy = 1- NGH/H %


BioPhase 8800 system with NF: Consistency and repeatability

Traditional UV detection: IgG control



NF detection IgG Control

NEW!



Comparable profiles observed with the same level of repeatability


BioPhase 8800 system with NF: Consistency and repeatability

Intra-Capillary	(LC+HC)%- UV		(LC+HC)%- NFD		HC/(HGH_HC)% -UV		HC/(HGH_HC)% -NFD	
	Average	% RSD	Average	% RSD	Average	% RSD	Average	% RSD
A	92.55	0.07	91.63	0.04	89.86	0.11	89.77	0.04
B	92.50	0.05	91.63	0.03	89.79	0.07	89.76	0.02
C	92.50	0.06	91.63	0.01	89.80	0.07	89.75	0.02
D	92.52	0.04	91.63	0.02	89.80	0.06	89.76	0.02
E	92.56	0.06	91.62	0.01	89.84	0.10	89.74	0.02
F	92.51	0.04	91.62	0.03	89.79	0.06	89.75	0.03
G	92.53	0.07	91.61	0.03	89.84	0.11	89.73	0.03
H	92.52	0.07	91.64	0.05	89.81	0.11	89.76	0.06
<i>Inter-capillary</i> <i>N = 48</i>	92.52	0.06	91.63	0.03	89.82	0.09	89.75	0.03

Comparable sample purity with the same level of repeatability


BioPhase 8800 system with NF: Emicizumab analyzed with CE-SDS

Non-reduced



LMW	Main	HMW
3.02%	96.85%	0.13%

Reduced



Glycosylation site occupancy: 98.01%



Sample: 0.5 mg/mL Emicizumab, Injection: 20 s @ -5 kV



AAV capsid protein analysis




AAV8 capsid protein analysis on BioPhase 8800 system



- Comparable capsid protein profiles across three different detections.
- Comparable repeatability across capillary, across injections was observed.

AAV8 capsid protein analysis on BioPhase 8800 system: dynamic range



The linearity plot of the CPA of VP3 versus the titer of an AAV8-CMP-GFP sample using the CE-SDS with UV 220nm, NF, LIF detection on the BioPhase 8800 system

	LDR [GC/mL]	LOQ [GC/mL]	LOD [GC/mL]	R ²
UV	4.7E+11 - 3.6E+13	9.4E+11	4.69E+11	0.999
NF	1.2 E+11 - 3.6E+13	4.7E+11	1.17E+11	0.999
LIF	9.2E+08 - 1.5E+13	9.16E+08	3.05E+08	0.999



AAV8 capsid protein analysis on BioPhase 8800 system: capsid ratio

Capsid Ratios

Experimental value			
	VP1	VP2	VP3
LIF	1.0	1.2	4.5
NFD	1.0	1.5	7.8
UV	1.0	1.4	7.5
Amino acid content			
Lys	32	22	17
Trp	15	12	12
aa	737	600	534
Normalized by amino acid content			
LIF	1.0	1.7	8.5
NFD	1.0	1.9	9.8
UV	1.0	1.7	10.4
Expected	1.0	1.0	10

VP1

AADGYLPDWLEDNLSEGIREWWALKPGAPKPKANQDQKQDDGRLVLPGYKYLGPFGDLKGEPVNAADA
AALEHDKAYDQQLQAGDNPYLRYNHADADEFQERLQEDTSFGGNLGRAVFQAKKRVLEPLGLVEEGAKT
APGKKRPVEPSPQRSPDSSSTGIGKKGQPARKRLNFQTGDSSESVPDPQLGEPRAAPS6V6PNTM

VP2

AAGGGAPMADNNEGADGVGSSSSGNWHCDSTWLGDRVITTSTRTWALPTYNNHLYKQISNGTSGGATNDN
TYFGYSTPWGYFDFNRFHCHFSPRDWQRLINNNWGFPRKRLSFKLFNIQVKEVTQNEGTKTIANNLTSTI
QVFTDSEYQLPYVLGSAHQGCLPPFPADVFMIPQYGYLTLNNGSQAVGRSSFYCLEYFPSQMLRTGNNFQ
FTYTFEDVPPFHSSYAHSQSLDRLMNPLIDQYLYLSRTQTTGGTANTQTLGFSQGGPNTMANQAKNWLP
GPCYRQQRVSTTTGQNNNSNFAWTAGTKYHLNGRNSLANPGIAMATHKDDDEERFFPSNGILIFGKQNAAR
DNADYSQVMTLSEEEIKTTNPVATEEYGVADNLQQQNTAPQIGTVNSQGALPGMVWQNRDVYLQGPWWA


VP3

KIPHTDGNFHPSPLMGGFGLKHPPPQILIKNTPVPADPPTTFNQSKLNSFITQYSTGQVSVEIEWELQKEN
SKRWNPFIQYTSNYUKSTSVDFAVNTEGVYSEPRPIGTRYLTRNL


Capsid ratio using NF is similar to the values using UV,
and closer to theoretical values

AAV9 capsid protein analysis on BioPhase 8800 system

scAAV9-CMV-LacZ



scAAV9-CMV-GFP



scAAV9-CMV-LacZ capsid ratio

	VP1	VP2	VP3
LIF	1.0	1.0	4.5
NF	1.0	1.0	8.1
UV	1.0	1.0	8.2

scAAV9-CMV-GFP capsid ratio

	VP1	VP2	VP3
LIF	1.0	1.2	5.2
NF	1.0	1.5	8.0
UV	1.0	1.4	8.6

Comparable trends for capsid protein ratios were observed when AAV9 is analyzed

AAV Capsid protein analysis using three different detection method

Signagen_AAVs-CMV-6FP	VP3				VP2		VP1	
	MT [min]	RSD% of MT	CPR%	RSD% of CPR%	CPR%	RSD% of CPR%	CPR%	RSD% of CPR%
LIF-600 nm	20.94	0.6	65.8	0.3	17.1	0.7	13.4	1.1
NF-360 nm	21.28	0.7	74.7	0.3	13.30	1.2	8.80	2.9
UV-220 nm	20.90	0.5	74.5	0.8	13.00	3.5	9.10	5.4

Signagen_AAVs-CMV-lacZ	VP3				VP2		VP1	
	MT [min]	RSD% of MT	CPR%	RSD% of CPR%	CPR%	RSD% of CPR%	CPR%	RSD% of CPR%
LIF-600 nm	20.50	0.4	70.30	1.2	14.00	1.90	13.40	1.1
NF-360 nm	20.89	0.8	81.34	0.6	9.48	1.27	9.18	4.0
UV-220 nm	20.38	0.6	80.85	1.1	10.06	5.76	9.09	6.6

Signagen_scAAVs-CMV-6FP	VP3				VP2		VP1	
	MT [min]	RSD% of MT	CPR%	RSD% of CPR%	CPR%	RSD% of CPR%	CPR%	RSD% of CPR%
LIF-600 nm	20.44	0.5	69.02	0.3	14.89	0.5	12.06	1.0
NF-360 nm	21.20	1.1	80.12	0.5	9.56	1.1	9.15	2.8
UV-220 nm	20.56	0.9	81.90	1.2	9.56	7.1	8.55	12.0


Comparable repeatability was observed between LIF, NF, and UV

A full-page photograph of a male scientist with a long white beard and safety glasses, wearing a white lab coat and blue gloves. He is using a pipette to transfer liquid into a clear bottle. The background shows a laboratory environment with various equipment and a poster on the wall.

cIEF analysis for proteins


BioPhase 8800 system with NF: cIEF of NIST mAb

UV detection: NIST mAb



NF detection: NIST mAb


NEW!



NF detection shows similar profiles for all charge variants with improved S/N ratios


BioPhase 8800 system with NF: cIEF of USP IgG

UV detection: USP IgG





NF detection: USP IgG

NEW!



NF detection shows similar profiles for all charge variants with improved S/N ratios

BioPhase 8800 system with NF: repeatability



Intra and Inter-capillary reproducibility. The inset shows the zoomed-out electropherograms.

BioPhase 8800 system with NF: Consistency and repeatability


Intra-Capillary	Acidic %		Main %		Basic %		Main pI Values	
	Average	% RSD	Average	% RSD	Average	% RSD	Average	% RSD
A	15.72	0.4	63.81	0.5	20.46	1.18	7.66	0.1
B	16.07	1.2	61.80	0.6	22.13	0.95	7.66	0.1
C	16.25	0.6	61.25	0.5	22.50	0.91	7.66	0.2
D	16.30	0.6	61.36	0.7	22.35	1.47	7.66	0.1
E	15.94	0.8	62.43	0.8	21.63	1.82	7.66	0.1
F	15.93	0.6	62.75	0.8	21.32	1.94	7.65	0.2
G	15.75	0.4	63.03	0.5	21.22	1.27	7.66	0.1
H	15.46	1.1	64.46	0.7	20.08	1.49	7.66	0.1
<i>Inter-capillary 48 injections</i>	<i>15.93</i>	<i>1.8</i>	<i>62.61</i>	<i>1.8</i>	<i>21.46</i>	<i>2.0</i>	<i>7.66</i>	<i>0.1</i>

Comparable sample purity with the same level of repeatability

cIEF Profile comparison between control and stressed sample

UV detection

Baseline drifts make integration on the basic side more subjective/uncertain




cIEF Profile comparison between control and stressed sample

NEW!

NF detection

More confident to integrate the basic side and reflect the changes

		Basic %	Main%	Acidic%
Control	UV	3.39	78.20	18.12
	NFD	3.39	78.55	18.06
Stressed	UV	6.28	71.84	21.81
	NFD	8.84	69.57	21.55




The data were obtained using a prototype BioPhase 8800 system with NFD



A full-page background image of a male scientist with a long white beard and safety glasses, wearing a white lab coat and blue gloves. He is using a pipette to transfer liquid into a clear bottle. The background shows a laboratory environment with various equipment and a poster on the wall.


CZE rapid charge variant analysis

CZE analysis of Trastuzumab on the BioPhase 8800 system and the PA 800 Plus system.



	S/N ratio [RMS] for the main peak		
	Trastuzumab		
Detection	NFD on BP	UV on BP	UV on PA
# of injections	42	42	6
Average	10658	8856	5533

CZE analysis of Trastuzumab on the BioPhase 8800 system with NF: repeatability



CZE analysis of Trastuzumab on the BioPhase 8800 system with NF: repeatability

	% Total basic variants			% Total main peak			% Total acidic variants		
Detection	NFD on BP	UV on BP	UV on PR 800	NFD on BP	UV on BP	UV on PR 800	NFD on BP	UV on BP	UV on PR800
* of injections	42	42	6	42	42	6	42	42	6
Average	7.5	8.3	9.2	72.9	71.3	71.2	19.6	20.4	19.6
STDEV	0.2	0.3	0.1	0.3	0.6	0.5	0.3	0.5	0.4
%RSD	2.5	3.1	1.6	0.4	0.8	0.7	1.7	2.3	2.2

Summary- NF detection for all protein workflows on BioPhase 8800 system

Native: signal relies on the naturally existing aromatic amino acid in the protein, eliminating the protein labeling process

Fluorescence: no direct transmission/absorbance interference from background media in the capillary, providing a clean background

BioPhase 8800 system: 8-capillary parallel processing and easy-to-integrate data, reducing the time to answer with high confidence

BioPhase 8800 system with NF : demonstrated to be suitable for all existing protein-related workflows, providing comparable results, faster and easier data processing.

Place holder for all the Sciex-related materials in CE Pharm.

Monday, Sep 8th

10:00-11:15 Monday, 8 September, 2025, Salon A-C

Session I - Application and Advancements in CE Innovations

Advancement of Capillary Gel Electrophoresis Through the Use of Native Fluorescence Detection

Jessica Taylor, *Teva Pharmaceutical Industries Ltd.*

11:45-12:45 Monday, 8 September, 2025, Foyer A-C

Lunch and Learn Technical Seminar Presented by SCIEX

Presentation type: LS - Livestreamed

Session Speakers:

Advancing Capillary Electrophoresis Sensitivity and Repeatability with Native Fluorescence Detection on the BioPhase 8800 System

Fang Wang, *SCIEX*

Label-Free Native Fluorescence for High-Sensitivity Protein Analysis using a Multi-Capillary Electrophoresis System

Brian Wei, *Sanofi*

Tuesday and Wednesday, Sep 9-10th

Sciex Booth #8

11:45-12:45 Tuesday, 9 September, 2025, Foyer A-C

44

Enhancing protein purity analysis using capillary electrophoresis-sodium dodecyl sulfate (CE-SDS) with native fluorescence detection (NFD) on the BioPhase 8800 system

[Tingting Li](#)¹, Jane Luo¹, Marcia Santos¹, Sahana Mollah²

¹ SCIEX, Brea, California, USA. ² SCIEX, Redwood city, California, USA

45

Assessing method transfer and equivalency between the PA 800 Plus system and the BioPhase 8800 system

Marcia Santos¹, Zoe Zhang², [Sahana Mollah](#)²

¹ SCIEX, Brea, California, USA. ² SCIEX, Redwood city, California, USA

46

Characterization workflows for CRISPR-Cas9 gene editing products, plasmids, and mRNA drugs on a single capillary electrophoresis (CE) platform

[Jane Luo](#)¹, Tingting Li¹, Kerstin Pohl², Sahana Mollah³

¹ SCIEX, Brea, California, USA. ² SCIEX, Marlborough, Massachusetts, USA. ³ SCIEX, Redwood City, California, USA

47

Advancing biopharmaceutical molecule charge heterogeneity analysis with native fluorescence detection on a multi-capillary electrophoresis platform

[Jane Luo](#)¹, Peter Holper², Marcia Santos¹, Tingting Li¹, Fang Wang³, kerstin Pohl³, Sahana Mollah²

¹ SCIEX, Brea, California, USA. ² SCIEX, Redwood City, California, USA. ³ SCIEX, Marlborough, Massachusetts, USA

67

Advancing Purity Analysis of Protein Therapeutics: Evaluation of the Next-Generation BioPhase 8800 System with Native Fluorescence Detection

Valdehi Bhagat, Maroly Gonzalez-Perez, [Douglas Johnson](#)

AstraZeneca, Gaithersburg, USA

49

Further Exploration of Native Fluorescence Detection in Capillary Electrophoresis

[Matthew Myers](#), Jessica Taylor

Teva Pharmaceuticals, West Chester, Pennsylvania, USA

34



Trademarks/ Licensing

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures. Waters Empower™ software is being used under license.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries [see www.sciex.com/trademarks].

© 2025 DH Tech. Dev. Pte. Ltd. MKT-35399-R



The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product[s] not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2025 DH Tech. Dev. Pte. Ltd. MKT-35399-A

