



TreeFrog
therapeutics

Cell therapy for all



Founded in 2018

Co-founded by biophysicist Kévin Alessandri & stem cell biologist Maxime Feyeux.
CEO Frédéric Desdouits appointed in 2021

100+ froggies

Blending biophysicists, stem cell biologists & bioproduction engineers.
HQ in Bordeaux + hubs in Boston and Kobe



\$82M raised

2021 - \$75M Series B; 2019 - \$7M Series A

bpifrance

LGP

LEONARD GREEN
& PARTNERS



Bristol Myers Squibb™

Building a global footprint in cell therapy

Two sites in Bordeaux, France



Clean rooms & analytics



Chip manufacturing & administration

Two international hubs to drive partnership opportunities in key markets



LabCentral 238, Kendall Square, Cambridge



Creative Lab for Innovation in Kobe (CLIK)

Unlocking iPSC-derived cell therapies




#1 Scale-up

Mass-produce billions of cells at low cost to serve thousands to millions of patients

#2 Quality

Secure cell quality at scale to ensure the safety of iPSC-derived cell products



#3 Clinic & Access

Improve access to treatments by lowering costs, providing new graft format with high efficacy

#1 Scale-up




Expected Benefit from iPSCs - Scaling up

Autologous	Allogenic	iPSCs
1 patient = 1 batch Patients own cells	1 donor or X donors = 1 batch = 1 to 50 - 100 patients	Banking strategy: 1 MCB = hundreds of vials, and by establishing WCB, you can in theory treat an unlimited number of patients.
BUT	Off the shelf availability	Off the shelf availability
Patients own cells	Material available for testing and validation	Off the shelf availability
Limited material available for testing and validation	BUT	Material available for testing and validation
No treatment if manufacturing failure	Batch to batch comparability may be challenging	BUT
Delay in administration	Donor/Patient compatibility*	Donor/Patient compatibility*
	*hypo-immune cells approach feasible	hPSCs specific challenges (see after)
		*hypo-immune cells approach feasible


C-Stem™ addresses bottlenecks for human PSC growth

outperforming 2D and 3D approaches while maintaining *in vivo*-like cell quality & functionality


Current approaches



2D COLONIES



AGGREGATES IN BIOREACTORS



C-STEM™ IN BIOREACTORS

	2D COLONIES	AGGREGATES IN BIOREACTORS	C-STEM™ IN BIOREACTORS
Yield (per week)	30- to 100-fold	10- to 40-fold	100- to 300-fold
Genomic risk	High Overcrowding and high mortality increasing selective pressure	Critical Overcrowding and high cell death (necrotic cells and impeller-induced cell death) driving selective pressure	Contained In vivo-like growth with low mortality; low selective pressure for oncogenic mutations
Scalability	Scale-out only Labor-intensive, batch-to-batch variability; high QC costs	Scale-up limited to 3L Due to impeller-induced shear stress damaging cells	No limitation in bioreactor size Cells protected by the capsule; cell culture parameters preserved across scales

High-throughput encapsulation

>1000 capsules / second



GMP encapsulation device

Automated, closed, single-use environment





A black and white micrograph showing a stream of small, dark, spherical capsules being dispensed from a nozzle at the top. The capsules are falling in a vertical column, with some scattered around. The background is light gray.

High-throughput encapsulation

>1000 capsules / second

64 hours




A micrograph of a single, teardrop-shaped capsule. Inside the capsule, there is a dense, circular cluster of cells. The capsule is surrounded by a thin layer of fluid.

in vivo-like
iPSC growth

with virtually no cell death

Alginate capsule

Highly reproducible exponential iPSC expansion in 10L bioreactors (2021)



>98%

High viability
(dead cells not rinsed away as in 2D)

>97%

Robust pluripotency across batches
95.2% OCT4 / 99.5% SOX2 / 98.3% NANOG

See Cohen *et al.*, Biomaterials, 2023

**Unprecedented iPSC amplification
in 10L bioreactors***

276x

in 6.59 days

1 batch, 1 week = 15Bn hiPSCs

*Best performance previously reported in literature: 37-fold in 6 days in 10L bioreactor (Huang *et al.* 2020), with significant drop in stemness

CONFIDENTIAL PRESENTATION

Cell Therapy for ALL

With the scale, cost reduction can be obtained by:



- Savings on materials and consumables
- Batch size allowing to treat hundreds of patients
- Automation of production
- Access to treatments despite limited starting material



#2 Quality

Benchmarking Cell Culture Systems for iPSCs

Continuous culture over 28 days



GENOMIC BENCHMARK


2D vs Aggregates vs C-Stem™

C-Stem: fastest growth, with reduced passages over 28 days




hiPSC quality : matching standard systems


(in which dead cells are rinsed away)




*Dead cells lost with media change



Controlling pre-oncogenic 20q11 mutation throughout large-scale expansion



Copy number values



Starting material: hiPSC cell line containing undetectable levels* of cells with a 20q11. 21 duplication. Head-to-head comparison of 2D, aggregates and C-Stem™ formats. While providing the best amplification performance, C-Stem™ is capable of controlling copy number variations throughout prolonged culture.

*below the limit of detection (LOD)


CONFIDENTIAL PRESENTATION

FOCUS ON: iPSC Cell Lines Requirements

- Genetic Stability is an ongoing issue in the pluripotent stem cell field
- Not All Cell Lines are of the same quality
- There is currently no standard approach on genes to be monitored for cell line selection



How good are iPSC Lines ?



Study of 839 hPSC lines, 285 failed testing (34%)

61% - Genetic Stability (G-banding)

- Half of which were recurrent abnormalities

25% - Viability/Usability

- No attachment, expansion or excessive differentiation

10% - Identity (Short Tandem Repeats)

- Identity mismatch, sex mismatch or mixed cell population

4% - Sterility

- Non-sterile or mycoplasma positive

D. Felkner., J. Brehm, E. McIntire, S. Minter, A. Paguirigan, K. Remondini, S. Taapken & T.E. Ludwig. "Human pluripotent stem cell quality: A scientific wake-up call." Poster presented at the ISSCR Annual Meeting; June 2019; Los Angeles, California.

What we know:


- **Karyotypic abnormalities are acquired throughout culture in a spontaneous and nonrandom manner (K 1, 12, 17, 20, and X)**
- **It provides competitive advantage to iPSCs through different mechanisms** (attachment after seeding, faster cell cycle times, reduced death, reduced differentiation...)
- **They are similar to karyotypic abnormalities observed in human cancers.**

Are iPSC Cell Line good enough for clinic?

Consensus = no go	In discussion	Unknown
Cell Lines with complete chromosome duplications or deletions cannot be used in clinic COSMIC list, PMDA List (2013) TP53, BCOR	What to do with SNV, CNV? ISSCR recommendations on Common recurrent genetic changes in hPSCs, some of their phenotypic consequences and suitable methods for their detection	What if below the level of detection? CNV Thresholds? Number? Size? Kinetics? Acceptability if no observation in NC and tumorigenicity studies? How does RBA support Clinical use

What's Next? if we want iPSCs' potential to become become a reality

- Industry:
 - Shared responsibility to look at the issue:
 - Cell Line suppliers
 - Drug Developers
 - Support the establishment of standards by generating data
 - Be transparent to support RBA if applicable.
- Regulators / Pharmacopeia:
 - Need a stakeholder workshop / White paper to pave the way for Guidance on iPSC Cell line quality requirements, choice is made very early in development (go/nogo)
 - Standards on oncogenic panel testing will harmonize practices, improve safety and allow for easier identification of the suitable materials from suppliers.




#3 Impact on Clinic & Access

C-Stem™: unlocking cell therapy potential




**hiPSC-derived
microtissues**



Neural
Lead program (Parkinson's)


in house



Liver
9 months from concept to IP

in house


partnering



Cardiac
6 months from concept to IP


in house

partnering



RPE
3 months from concept to IP

in house




CD34+
6 months to launch *in vivo* PoC


partnering




**Donor-derived
microtissues**



T cells
3 months from concept to IP




NK cells
On-going PoC




MSCs
6 months for PoC study

> 100 hundred trials based on pluripotent stem cells yet, no significant clinical progress beyond Phase I/II

Number of clinical trials based on pluripotent stem cells launched per year



Number of clinical trials based on pluripotent stem cells per disease area






Parkinson's disease cell therapy

C-Stem™ -based GMP Manufacturing Process & QC

Bringing automation, 800 doses in 1L bioreactor



42 days

In-process QC

- Viability
- Identity
- Morphology

Batch release QC

- Quantity/Viability
- Purity
- Impurities
- Potency (dopa release)
- Karyotyping


Additional QC

- Genomic integrity (CNV)
- Single-cell characterization

Neural microtissue engrafts, dopaminergic neuron projections observed

The implanted human microtissues survived, thrived in the lesioned hemisphere of the rat brain post implantation

TH positive dopaminergic neurons cell bodies detected in the graft area



Cryopreserved microtissues transplanted (high dose) within the dorsal striatum after a 5 months $\pm$ 10 days of *in vivo* experiment period
(A) Overview of the whole grafted area and (B) observations of the TH+ cells projection within the host tissue.

Full motor function recovery: 8 weeks



Hemiparkinsonian pre-graft (rat #13)

Hemiparkinsonian pre-graft (rat #13)



8 weeks after C-Stem™ neuron implantation (rat #13)

C-stem™ platform : unlock iPSC potential



Fulfilling
the iPSC promise
From 1 good edit
to billions of cells



Meet CMC expectations
for allogeneic products Scale,
automation,
quality



Create unique
products for various indications with
unique phenotypic specs
unique 3D format



TreeFrog

therapeutics

Cell therapy for all



@treefrogtx TreeFrog



Stem Cell Jungle

www.treefrog.fr

francois.gianelli@treefrog.fr
michael.lanero@treefrog.fr