

WCBP 2023 Waters Technical Seminar

Nick Pittman Marketing Manager, Global Biopharmaceutical Business 26th January



Protein Digestion Quick, Clean, Complete

Challenges with Protein Digestion





Protein digestions are plagued by *long digestion times and performance trade-offs*

Productivity and data analysis can be hampered by messy maps with various non-tryptic peaks

What if you had reproducible, efficient digestions that deliver high sequence coverage and low baseline noise?

Quick, Clean, Complete





MS-Grade Trypsin

Coming February 2023...

QUICK | Tired of long sample prep times?



High enzyme ratios (1:5) are faster, but typically lead to high autolysis. RapiZyme trypsin's unique autolysis resistance and low missed cleavage unlocks *efficient 30-minute digestions* vs. standard 3-hour digestion (using 1:20)

Remicade[™] digestion, 1:5 Enzyme: Protein Ratio, 30 min.

Trypsin Autolysis (% of TIC)

Missed Cleavage (% of TIC)



Digestion details: 37 °C, pH 7.5. N=6 per enzyme, 2 batches of each enzyme, 3 digestion replicates per batch.

CLEAN | Looking for noise-free maps for quick and easy data analysis?

RapiZyme Trypsin <u>% Impurity 1</u> 12.3% 30 37.5 2.5 7.5 10 12.5 15 17.5 20 22.5 25 27.5 32.5 35 40 42.5 45 5 3.8% Leading Competitor High Background Noise Autolysis peaks Leading Competitor RapiZyme 37.5 2.5 7.5 10 12.5 15 17.5 20 22.5 25 27.5 30 32.5 35 40 42.5 45 Retention time [min]

Critical decisions are made distinguishing the smallest changes in a peptide map. Proceed quickly and confidently with **exquisite baselines** free from noise

1. % total TIC area associated with autolysis, unmatched, non-specific, and missed cleavage peak area. N=6 digestion replicates, 2 batches of each enzyme, 3 digestion replicates per batch. Chromatogram is of 1:5 Enzyme:Protein ratio digestions of Remicade, 30 minutes, 37 °C, pH 7.5.

% 75·

50·

25

100-

%

75·

50·

25

Waters™

COMPLETE | Need to ensure full identification and characterization?

<u>% Sequence Coverage</u>



<u>1:5, 30 Minutes Digestions</u>: Remicade[™] digestion. Average of n=6 per enzyme, 2 batches for each enzyme, 3 digestion replicates per batch, 37 °C, pH 7.5 <u>Standard Digestion</u>: 1:20 enzyme:protein ratio, 3 hours, 37 °C, pH 7.5, leading competitor enzyme, 1 batch, 1 digestion replicate

Waters™



How Can You Ensure Maximum Confidence in Your Results?

What Happens When You DON'T Have Confidence in Your Results? Waters™

You inject a sample But see this.... and expect this.... 776.66 > 94.9 (GEM91 B) 7.29e4 MRM of 4 Channels ES-776.66 > 94.9 (GEM91 B) 100 481 Injection 3 1.00 2.00 3.00 4.00 1.00 1.20 1.40 0.20 0.40 0.60 0 80 1.60 1.80 2.00 2.20 2.40 2.60 2.80 3.00 3.20 3 40 3 60 3.80 776.66 > 94.9 (GEM91 B) 4.26e4 MRM of 4 Channels ES-776.66 > 94.9 (GEM91 B) 100 Injection 2 1.00 2.00 3.00 4.00 0.20 0.40 0.60 0.80 1.00 1.20 1.40 1.60 1.80 2.00 3.00 2.20 2.40 2.60 2.80 3.20 3.40 3.60 3.80 776.66 > 94.9 (GEM91 B) 2.21e4 MRM of 4 Channels ES-776 66 > 94 9 (GEM91 B) 100-325 Injection 1 1.00 2.00 3.00 4.00 0.20 0.40 0.60 0.80 1.00 1.20 1.40 1.60 1.80 2.00 2.20 2.40 2.60 2.80 3.00 3.20 3.40 3.60 3.80

Where is it?

100

100

100

a

Molecules Can Bind to Surfaces



Non-Specific Binding (NSB), Non-Specific Adsorption (NSA) can be a problem

- Any binding or adsorption that was NOT intended
- Molecules tend to adhere to any exposed surfaces
- Any chemical interaction can be the source of binding, but especially:
 - Polarity-based interactions (e.g., hydrophobic attraction)
 - Ionic interactions (e.g., coulombic attraction)

To suppress binding losses of known interactions...

- Avoid interaction between the surface and analytes
- Maintain the environment so that interactions are not strong



Some LC Separations are More Challenging Than Others

Biomolecular Analytes Prone to Metal Interactions



Waters[™]

A Better Solution from Waters R&D

Waters[™]

HIGH PERFORMANCE SURFACES



COLUMNS AND SYSTEMS

with MAXPEAK







More Information and Additional Educational Assets

Waters™



SEC Premier Classes of therapeutic proteins separated on ACQUITY™ Premier Protein SEC 250 Å, 1.7 µm Column ADC mAb Size Variant Std Kadcyla Column: Dimension: Flow rate: Temp.: Wavelength: Mobile phase: Ado-trastuzumab emtansir Injection Load: 12 µg 4.6 x 300 mm 0.3 mL/min 35 °C Instrument: 8.00 9.00 tc.00 7.00 min 8.00 9.00 6.00 0.0018 Fusion Fc Fusion mAb Orencia (Abatacept) Injection Load: 12.5 µg laG1 NistmAb RM 8671 Injection Load: 7 µg Protein 0.0012 7.00 min 8.00 9.00 9.00 10.00 7.00 min 8.00 6.02 6.03 lgG1 Chimer 0.0009 lgG1 0.0036 IgG2 Rituximab Injection Load: 7 µg Panitumumab Injection Load: 12.5 µg 7.00 min 8.00 6.00 7.00 min 8.00 9.00





Peptide Premier

Glycan Premier OUT-OF-THE-BOX PERFORMANCE Conventional MaxPeak Premi Technology Passivated with Fetuin Passivated with Fetuir 1500000.12 1000000.0 1000000.0 FA2BG2S 500000.0 500000 lni #4 Ini #4 2500000. 2000000.1 1500000 12 1000000.0 1000000.08 500000.0 500000 lnj #1 lnj #1

EADBOOS



Performance from Injection #1 with RapiFluor-MS™ Glycan Performance Test Standard



bioseparation



waters.com/MaxPeakColumns waters.com/BioAdvisor

©2023 Waters Corporation



Waters @ WCBP 2023



New Literature

MAM Application Notebook

Collation of Waters MAM material for BioAccord & QDa <section-header>

Multi-Attribute Methods for Biopharmaceutical Analysis



MAM eBook Interviews with subject matter experts







Chromatographic methods for the analysis of RNA therapeutic compounds; Matt Lauber, Joe Fredette

Evolving Biotherapeutic Methods from Characterization to Commercialization; *Robert Birdsall,* **Scott Berger**

mponents can be identified with further data analysi





REGENERON

Implementation of a Qualified HRMS System in the QC Environment

Yuwei (Serena) Wu