

An abstract graphic featuring a complex network of interconnected nodes and lines, resembling a molecular structure or a data network. The nodes are represented by small circles in various shades of blue and grey, connected by thin, light blue lines. The background is a gradient of blue, with a darker horizontal band across the middle where the text is located.

Waters lunchtime seminar, WCBP 2025

Partnering with Waters for Biopharmaceuticals analyses: Innovations supporting product development and release

Nick Pittman

*Marketing Manager,
Global Biopharma Business*

Life saving therapies

Waters™

Breast cancer patient survival after 4 years¹

85%

with
immunotherapy

67%

only
chemotherapy



Crohn's disease patient remission after 1 year²

36%

with
immunotherapy

12%

with
placebo



Multiple Sclerosis patient free of lesion after 1 year³

93%

with
immunotherapy

19%

with
placebo




Sources: ¹Rodmon EH et al. N Engl J Med. 2005 ; Slamon DJ et al., N Engl J Med. 2001 ; ²www.humira.com/crohns/about-humira/results-with-humira ; ³Yamout et al., J. Imm. Res. 2018

Developing biopharmaceuticals is a tough business


COST

Average R&D cost to develop a therapy (Discovery to launch)




TIME

Enrolment duration for industry sponsored clinical trials (months)




RISK

Composite success rates (from IND to approval)



VALUE

Average peak product sales (Per pipeline asset)



More efficient working practices are needed to mitigate these headwinds

How to foster access to life saving therapies to all those in need?

Waters™

Faster and cost-effective development and production of biotherapeutics

3 key industry drivers

COST



1 Confidence

Deliver the right results, first time.
Meet requirements from Discovery,
to Development & QA/QC

PRODUCTIVITY



2 Consistency

Continued technology evolution
for the most effective & reliable
separations & detection

COMPLIANCE



3 Connectivity

Address critical regulatory needs
with seamless, streamlined LCMS
& Light Scattering workflows

How to foster access to life saving therapies to all those in need?

Waters™

Faster and cost-effective development and production of biotherapeutics

3 key industry drivers

COST



1 Confidence

Deliver the right results, first time.
Meet requirements from Discovery,
to Development & QA/QC

PRODUCTIVITY



2 Consistency

Continued technology evolution
for the most effective & reliable
separations & detection

COMPLIANCE



3 Connectivity

Address critical regulatory needs
with seamless, streamlined LCMS
& Light Scattering workflows

1 Confidence

Instrument solutions deliver results at point-of-need
and meet specific requirements across the product lifecycle

Waters™

Discovery

Development

Production / QC

Non-regulated

Regulated / GxP

Little product knowledge: **Characterization**

- Leading edge analytical instruments
- Expert level

Good product knowledge: **Product and process monitoring / QC**

Accessible, rugged and robust -
Routine and throughput -



Andrew[™]
the pipetting robot



Core technologies

1
Confidence

2 New instrument solutions deliver results at point-of-need Waters™ and meet specific requirements across the product lifecycle

Discovery

Development

Production / QC

Non-regulated

Regulated / GxP



Xevo™ | MRT



Acuity™ QDa II

1 Confidence

2 New instrument solutions deliver results at point-of-need Waters™ and meet specific requirements across the product lifecycle

Discovery

Development

Production / QC

Non-regulated

Regulated / GxP



Xevo™ | MRT



Acuity™ QDa II

1 Confidence

Xevo MRT: High sensitivity multi reflecting time-of-flight MS

Increased assurance of reproducibility and accurate assignments



Exceptional
performance & speed
without compromise



Resolution
100K (FWHM)
@ m/z 956

Speed
50 Hz MS
100 Hz MS/MS

Sensitivity
MS/MS
x10 increase*

Mass Accuracy
<500 ppb
RMS

TIC



Rep.	Cmpnt	Obs. mass (Da)	Mass error (ppm)	Cmpnt (Rel %)
1	Fd'	25688.8476	-0.58	76.78
2	Fd'	25688.8616	-0.04	77.31
3	Fd'	25688.8503	-0.48	77.69
4	Fd'	25688.8575	-0.20	77.27
5	Fd'	25688.8569	-0.22	75.10
6	Fd'	25688.8538	-0.34	77.03
7	Fd'	25688.8545	-0.31	75.86
8	Fd'	25688.8584	-0.16	74.97
9	Fd'	25688.8711	0.33	77.26
10	Fd'	25688.8574	-0.20	77.57
Avg		25688.86	-0.22	76.69
RSD		0.00		1.31
STD.dev		0.01	0.25	1.01


Reliable insight and confident data in early development can inform project decisions when change costs are low

1 Confidence

Xevo MRT: peptide mapping data

Increased assurance of performance, reproducibility and assignments

Waters™



98% sequence coverage

Peptides assigned using max 2 ppm mass error

Linearity: 2000 ng - 0.02 ng on column


High sensitivity & linear MS response enables confident tracking of low-level CQA's across samples

1 Confidence

Xevo MRT: synthetic oligonucleotide data

Increased assurance of performance, reproducibility and assignments


Waters™



Results TIC TUV 260 MS No Calculation Results

Identity: **Pass** Purity: **Pass**

	Component	Observed mass (Da)	Expected mass (Da)	Mass error (ppm)	Observed UV RT (mins)	LC amount (%)	MS amount (%)	MSxLC (%)
2	7,770.6809	7,770.687	7,770.681	0.8	4.40	89.8	70.0	62.9
3	7770.6809 PS to PO	7,754.673	7,754.704	-3.9	4.40	89.8	5.6	5.1
4	7770.6809 CNET	7,823.703	7,823.707	-0.6	4.80	10.2	42.3	4.3
5	7770.6809 +TEA	7,871.768	7,871.801	-4.3	4.40	89.8	4.4	4.0
6	7770.6809 +TEA(2)	7,972.897	7,972.922	-3.1	4.40	89.8	1.3	1.2
7	7770.6809 +TEA,Oxidation M	7,887.771	7,887.796	-3.2	4.40	89.8	0.7	0.6
8	7770.6809 +TEA,-CH3	7,856.786	7,856.778	1	4.40	89.8	0.6	0.5
9	7770.6809 +Na,+CH3	7,807.691	7,807.686	0.5	4.80	10.2	2.7	0.3
10	7770.6809 +TEA,CNET	7,924.805	7,924.828	-3	4.80	10.2	2.7	0.3
11	7770.6809 -CH3	7,755.668	7,755.657	1.3	4.80	10.2	2.4	0.2



INTACT **Mass**

Identity confirmed within a mass error of 1 ppm

Efficient MS^E fragmentation confirms sequence

Automatically confirm FLP sequence and confirm the presence and identity of impurities

1 Confidence

2 New instrument solutions deliver results at point-of-need Waters™ and meet specific requirements across the product lifecycle

Discovery

Development

Production / QC

Non-regulated

Regulated / GxP



Xevo™ | MRT




Acuity QDa II

Evolution of the ACQUITY QDa™ II


Smart Innovation for a solution widely deployed across development, manufacturing, and quality control

Extend dynamic range below UV baseline
Follow molecules with no chromophores
Existing QDa methods directly transfer


- ✓ Extended mass range to 1500 m/z
- ✓ Improved repeatability and robustness
- ✓ Up to 60% Reduction in Service Time
- ✓ Up to 70% Energy Savings vs. SQ systems



Wide Mass Range



+35% increase in MS response




All CDR peptides detected for successful ID test. Extended mass range enables increased MS response

1 Confidence

Application example: monoclonal antibodies

Reproducible and consistent response at low level is required for impurity monitoring

Waters™



Sensitivity and consistent peak area shows suitability for quantitative assays that demand precision and accuracy

How to foster access to life saving therapies to all those in need?

Waters™

Faster and cost-effective development and production of biotherapeutics

3 key industry drivers

COST



1 Confidence

Deliver the right results, first time.
Meet requirements from Discovery,
to Development & QA/QC

PRODUCTIVITY



2 Consistency

Continued technology evolution
for the most effective & reliable
separations & detection

COMPLIANCE



3 Connectivity

Address critical regulatory needs
with seamless, streamlined LCMS
& Light Scattering workflows

2 Consistency

Continued technology evolution

MaxPeak Premier family expands: Alliance™ iS Bio HPLC System

Waters™

Eliminate unpredictability of analyte losses due to metal interactions

Ultimate performance for more consistency & confidence in analytical data

- ✓ Increase sensitivity & sharpen peak shapes
- ✓ Minimize variability of results
- ✓ Reduce passivation & conditioning times



LC technology




alliance[™]is | bio

Useability and transferability

- Reduce human errors up to 40%
- Up to 3X faster method migration and operator training

Improved, robust performance

- Lower carryover
- Increased injection precision
- Inert for Bioseparations



A routine system that introduces risk mitigation for enhanced productivity

2 Consistency

Continued technology evolution

Transforming SEC Analyses for Peptides & Proteins

Waters™



Column
hardware




US and International Regulatory Agencies require submission of precise and accurate size-variant separations of biotherapeutic peptides (e.g., GLP-1 receptor agonists) and small protein drugs

Why 125Å SEC MaxPeak Premier Columns?

Compared to USP Listed SEC Column for insulin analyses:

- **+15% Resolution** on HMWS vs Insulin Monomers
- **6x Faster** than Current HPLC-based Methods
- **7x Organic Solvent Savings**
- **Unmatched "Out of Box" Performance** for enhanced assurance for reliable data

1.7 and 2.5 μm particles for seamless scalability and method transfer



Consistent SEC separations help boost productivity & minimize risk on regulatory submissions & release testing

2 Consistency

Continued technology evolution

Advances for peptide and oligonucleotide purification

Waters™




Column
hardware



Why MaxPeak Premier purification columns?

- **Predict Scale-Up:** Transition seamlessly from scouting methods to bench scale preparative columns, ensuring consistent scale up results
- **Protect Sensitive Compounds:** **Inert** preparative columns provide extra assurance for recovering metal-sensitive compounds such as peptides and oligos from complex mixtures
- **Boost Isolation Efficiency and Productivity:** Better resolution and/or reduced preparative separation times with sub-5 um sorbent packing



Enhanced productivity through seamless, scalable purification for challenging therapeutic compounds

How to foster access to life saving therapies to all those in need?

Waters™

Faster and cost-effective development and production of biotherapeutics

3 key industry drivers

COST



1 Confidence

Deliver the right results, first time.
Meet requirements from Discovery,
to Development & QA/QC

PRODUCTIVITY



2 Consistency

Continued technology evolution
for the most effective & reliable
separations & detection

COMPLIANCE




3 Connectivity

Address critical regulatory needs
with seamless, streamlined LCMS
& Light Scattering workflows

3 Connectivity

Evolution of Wyatt technology into Waters portfolio

Bringing connectivity and compliance for biologic modalities



Adopt adjacent technology



The 2023 **acquisition of Wyatt Technology** has allowed Waters to build synergies to better support our customers

**Digital connectivity through HPLC
CONNECT with ASTRA™ software and
HPLC + UPLC instruments**



Keep up with complexity



Innovation of biopharmaceutical tools for new modality analyses

ZetaStar DLS/ELS/SLS Instrument + Arc HPLC

GTxResolve Column Chemistries, 1000A, 125A verified by MALS and **application development**

Future: Focus on regulatory needs



Demand for deeper characterization increasing with complexity

Migrate more analytical workflows into Empower to strengthen data integrity position

Streamlined analysis, data collaboration, enhanced data management


Deeply understand growing complexity of biologics with customers, delivering next generation connected, compliance ready solutions

3 Connectivity

Digested oligo mapping: enabling routine analysis

Simple workflow tools for confident sequence assignment

Waters™



*MS^e characterization with  CONFIRM **Sequence** or third-party data connectivity

Compliance – ready workflows for monitoring of RNA sequence through development and manufacture

Connectivity

New enzymes for flexible and optimized RNA digestion Waters™
sgRNA digestion patterns for common RNases

RNase T1¹

5' mG*mA*mU*G AUCUCUCAACUUUAAACGUUUAGAGCUAGAAAAUAGCAAGUUAAAAUAAAGGCUAGUCGCUUUAUCAACUUGAAAAAGUGGCACCGAGUCGGUGCmU*mU*mU*U3
mG*mA*mU*Gp UUUUAGp AAAUAGp UUAAAAUAGp UCCGp AAAAGp CmU*mU*mU*Up

RNase 4¹

5' mG*mA*mU*GAUCUCUCAACUUUAAACGUUUUAGAGCUAGAAAUAGCAAGUUAAAAUAGGCUAGUCCGUUAUCAACUUGAAAAAGUGGCACCGAGUCGGUGCmU*mU*mU*U3'
mG*mA*mU*GAUCUCUCAACUUU AAGCUUUUp AGAAAU AAAAU AGUCCGUU GAAAAAGU GCmU*mU*mU*U
 AACGUUUU AGAAAUp AAAAUp AGUCCGUUp GAAAAAGUp
 AACGUUUU AGCAAGUUp AAGGCU AUCAACUU GGCACCGAGUCGGUp
 AGAGCU AAGGCUp
 AGAGCUp

5' mG*mA*mU*GAUCUCUCACACUUUUAACGUUUUAGAGCUGAGAAAUAAGCAAGUUUAAAAUAAGGCUAGUCCGUUAUCAACUUGAAAAAGUGGCACCGAGUCGGUGCmU*mU*mU*U3'
 UCUCACAU UAACGU UAGAGC UAAGAA UAAGCCGU UGAAAAAGUGGCACCGAGUCGGUGC
 UCUCUCAC UUAACGU UAGCAAGU UAAGCCU UAGUCCGU UUGAAAAAGUGGCACCGAGUCGGUGC
 mG*mA*mU*GAUC UCAAC UAACGUU UAGAGCUGAGAA AAGU UAAAAUAAGGC UAUC UAUACAAC ACCGAGUCGGUGC
 UCAACU UUAAGAGC UAAGCC UAUAACU GAUC GGUGC
 UCUCAC UUUAGAGC UAGUCCGU UAUAACU mU*mU*mU*U
 UCAACUU
 AACU
 UUUAAAGU
 UUUAAAGUU
 UUUAAAGUUU
 UUUAAAGUUU

5' mG*mA*mU*GAUCUCUCAACUUUUAACGUUUUAGAGCUGAGAAUAGCAAGUUUAAAAUAAGGCUAGUCCGUUAUCAACUUGAAAAAGUGGCACCGAGUCGGUGCmU*mU*mU*U3'
 mG*mA*mU*GAUC AUC UUAAGAGC AGAAUAGCAAGUU UAAGGU GUUAUCAAC AUC GAGUC mU*mU*mU*U
 mG*mA*mU*GA UCUCAAC AGAGC AAAU UAGUCC UAUCAAC UGGCACCAGUC GGUGC
 UCAAC UAGAGC UAAAA UAUC AACUUGAAAAAGUGGC
 AACUUU UUAAGAGC UAGCAAGU UAGAAAA
 UUUUAA UAGAGC UAGAGC AAAA
 UUAACGUU UAGAGC UAGAGC
 AACGUU UAGAGC UAGAGC
 AACGUU UAGAGC UAGAGC
 UAGAAA UAGAAA

METMCUSATIVIN

Longer digestion products with greater sequence overlap and increased redundancy

Longer digestion products reduce ambiguity, allowing differentiation at MS1 level

How to foster access to life saving therapies to all those in need?

Faster and cost-effective development and production of biotherapeutics

Waters™

3 key industry drivers

COST



1 Confidence

Reduce costs with confident decision making

xevo™ | MRT Acuity QDaII

PRODUCTIVITY



2 Consistency

Reduce troubleshooting and investigations

MAXPEAK PREMIER alliance™ is bio

COMPLIANCE



3 Connectivity

Seamless workflows to meet regulatory needs


Waters™ | WYATT TECHNOLOGY | waters_connect

Waters | Wyatt Technologies at WCBP 2025: Come & see us at Booth #27

Lipid Nanoparticle Impurity Analysis: Leveraging Mass Spectrometry to Monitor Impurities and Aldehyde in Raw Materials and Formulated Drugs in Regulated Environments
Quintan Diaz, Robert E. Winkler
Waters Corporation, Milford, MA

CONNECTING SCIENCE
to Solve Problems, Fast


Waters™



Complementary Light Scattering Solutions for Protein Characterization


Bud Halim, Colette Quinn, Sophia Kenrick, Michelle Chen | Wyatt Technology, LLC

Waters™ | WYATT TECHNOLOGY



MEETING DEMANDS FOR IMPROVED MS PERFORMANCE TO SUPPORT DEEPER PRODUCT LC-MS CHARACTERIZATION WITH MULTI-REFLECTING TOP MS TECHNOLOGY


Jonathan Fair, Scott Berger, Nick Peterson, Ying Qiang Yu
Mettler Corporation, Waltham, MA, Waters Corporation, Milford, MA



TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

IMPROVING THE COVERAGE AND QUALITY OF RNA DIGEST OLIGONUCLEOTIDE MAPPING USING LC-MS, NOVEL DIGESTION ENZYMES, & AUTOMATED INFORMATION DATA PROCESSING

Cathy Conner, Bob Adams, Naveed Anjum, Chris Pinner, Matthew Goodwin, Jaka Bliznjak, Scott Berger, and Ying Qiang Yu
Mettler Corporation, Waltham, MA, Waters Corporation, Milford, MA



TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

Virtual Posters:

#51 Lipid Nanoparticle Impurity Analysis

#54 Complementary Light Scattering Solutions for Protein Characterization

#61 Deeper product understanding with novel MRT-TOF technology

#62 sgRNA / mRNA oligo digest mapping



Pre-conference webinar available On Demand:
Light Scattering-Based Solutions for Excipients and Advanced Formulations

Luncheon technical presentation:

Integrating LC, MALS, and
MS for Compliance and
Connectivity in GMP Labs



Serena Wu
Senior Analytical Specialist
Regeneron




Integrating LC, MALS, and MS for Compliance and Connectivity in GMP Labs

Serena (Yuwei) Wu

29 Jan 2025

REGENERON®



Agenda



- Transferring from Astra 6 and Astra 8: SEC-MALS for Reference Standard Characterization
 - Instrument/Software Comparison
 - Data Equivalency Study
- High Resolution Mass Spectrometer: routine release testing
 - Peptide Mapping ID Method
 - Alternative Format ID Method
 - MAM Capabilities

Agenda

- **Transferring from Astra 6 and Astra 8: SEC-MALS for Reference Standard Characterization**
 - **Instrument/Software Comparison**
 - **Data Equivalency Study**
- **High Resolution Mass Spectrometer: routine release testing**
 - Peptide Mapping ID Method
 - Alternative Format ID Method
 - MAM Capabilities

Background information on SEC-MALS

Molar mass of conjugated protein and modifier is measured based on data gathered from multiple detectors after separation



Size Exclusion Chromatography

UPLC or
HPLC:
UV Detector;
Concentration
Detector
+
SEC Column


Operated on Empower

Multi-Angle Light
Scattering Detector

Astra 6 (current)
Astra 8 (new)



RI Detector

Refractive Index; Concentration Detector



RS Characterization Assay: SEC-MALS

- **What is Reference Standard (RS)?**
 - It is defined as a batch of drug substance or formulated drug substance selected and qualified to ensure accurate, reliable and consistent analytical measurement.
 - “A well-defined home-made ruler to measure the quality of Regeneron products”
- **Why perform extended characterization on RS?**
 - FDA and ICH guidelines requires additional characterization procedures beyond routine testing for reference material.
 - “We need to prove that our ruler is well-defined before measuring everything against it”

Characterization Method	SEC-MALS
Quality attributes	Molar mass of conjugated protein and modifier
Current platform	Astra 6 with Empower FR2
 New platform	 Astra 8 with Empower FR4

Comparison between Astra 6 vs. Astra 8 platforms

Updated instrument platform brought enhanced data integrity control



miniDAWN
(TREOS)

Optilab (T-rEX)



miniDAWN
(NEON)

Optilab

Astra 6

1. System approaching end-of-life and lacks continued vendor support by the end of 2024.
2. Runs on Empower FR2 which is being decommissioned.
3. Refractive index of mobile phase needs to be entered manually.
4. User permission levels not separated.
5. Detector front panel could be accessed during data collection.

Astra 8

1. Continued vendor support under service plan.
2. Compatible with Empower FR4.
3. Refractive index model built in with instrument method.
4. Different permission level for analyst and method developer.
5. Detector front panel is now locked during data collection.

Qualification Requirements

1. System qualification

- Creation of operation and processing SOPs, data review WI
- Assist with data integrity findings

2. Method qualifications

- Transferring current CTPs to ASTRA 8
- Method quality of life updates


System Comparisons

Equivalency Study between Astra 6 and Astra 8




Astra 8 Accuracy

Molar mass from Astra 8 was compared against target molar mass value



Data equivalency between Astra 6 and Astra 8

Molar mass from Astra 8 and historical qualification data from Astra 6 were compared for %RSD



Section Summary

ASTRA 6

- End of life
- Data integrity
- Elevated privileges for analysts
- Certain processing parameters need updating

ASTRA 8


- Data integrity solutions
- Opportunity to update methods for best practices
- Opportunity to simplify procedures for analysts
- Proven comparability of molar mass data between systems

Agenda

- Transferring from Astra 6 and Astra 8: MALS for Reference Standard Characterization
 - Instrument/Software Comparison
 - Data Equivalency Study
- **High Resolution Mass Spectrometer: routine release testing**
 - **Peptide Mapping ID Method**
 - **Alternative Format ID Method**
 - **MAM Capabilities**

Implementing Mass Spec in QC Setting

Mass range:
50-7000 m/z for positive ions
50-5000 m/z for negative ions
Mass resolution: 10,000 FWHM




Mass range of up to 16,000 m/z
Mass resolution: 30,000 FWHM



- Protein Identification
- Protein MAM
- Impurities
- Amino Acids
- Glycans
- Reference Standard Characterization
- Oligonucleotides Identification
- Oligonucleotide MAM
- Virus-like particles
- Adeno-associated virus

Peptide Mapping Identification

- Per ICHQ6B, all identity tests of drug substances must be highly specific and should be based on unique aspects of its molecular structure and/or other specific properties
- Unique CDRs gives a unique UV profile used for identification
- New monoclonal antibodies with engineered modifications
 - 3 amino acid difference



Advantages of LC-MS for Peptide Map

LC-MS Peptide Map Method

Traditional UV Assay

- Long analysis times (2 hours/sample)
- Low throughput
- Ambiguity in data analysis

BioAccord:10 min method

- Short analysis times (10 mins/sample)
- High throughput
- Simplified yes/no data analysis output



Challenges in robustness and repeatability:

- Changes in peak resolution
- Peak co-elution or splitting

- Certain of peptide ID
- Changes to peak resolution and splitting have no impact

Demonstrating Specificity of the 10min peptide method

Testing scientific library of that's specific to mAb tryptic digest standard against a Regeneron antibody

mAb tryptic digest standard:
ID acceptance criteria **met**

Peptide Sequence	Chain	Expected mass (Da)	Identification Status
ALEWLADIWWDDK	HC	1660.8006	Identified
LTISK	HC	561.3606	Identified
VTNMDPADTATYYCAR	HC	1848.7891	Identified
DTLMISR	HC	835.4342	Identified
GFYPSDIAVEWESNGQPENNYK	HC	2544.1314	Identified
VTITCSASSR	LC	1081.5306	Identified
SGTASVVCLLNNFYPR	LC	1797.8952	Identified
VDNALQSGNSQESVTEQDSK	LC	2135.9687	Identified

Antibody 1: ID acceptance
criteria **not met**

Peptide Sequence	Chain	Expected mass (Da)	Identification Status
ALEWLADIWWDDK	HC	1660.8006	Identified
LTISK	HC	561.3606	Not observed
VTNMDPADTATYYCAR	HC	1848.7891	Not observed
DTLMISR	HC	835.4342	Identified
GFYPSDIAVEWESNGQPENNYK	HC	2544.1314	Identified
VTITCSASSR	LC	1081.5306	Not observed
SGTASVVCLLNNFYPR	LC	1797.8952	Identified
VDNALQSGNSQESVTEQDSK	LC	2135.9687	Identified

Summary

LC-MS Peptide Map Method



- Implements platform method
- Increases selectivity and specificity between nearly identical proteins
- Generates compliant-ready output
- Reduces ambiguity of data processing

- Automated sample preparation process greatly increases throughput
- Mitigates human error and reduce work-related injuries


Implementing MAM for Improved Efficiency

Challenge

- Conventional techniques require several assays for PQAs
- Some assays require extensive prep and have long run times to ensure resolution

Goal

- Develop methods with improved efficiency using LC-MS capable of detecting and quantifying several PQAs
 - Shorten run times
 - Eliminate repeat testing



MAM Capabilities by 10 min LC-MS Peptide Mapping

Product
identification by
coverage map




Create scientific
library for
peptides of
interests




Process by
accurate mass
screening

- FDA Information Request – Product oxidation study in FEP bag

Extracted ion chromatogram of Oxidized Met containing peptide





• Bioaccord Peptide MAM App



- Monitors PTMs across different samples and timepoints
- Calculates percent modification


Sample #	Condition
1-6	Oxi day: 0, 7, 14, 21, 28, 42
7-12	Dea day: 0, 7, 14, 21, 28, 42
13	Standard



Summary: Future of MAM


- More powerful detector gives potential for MAM
- Qualitatively and quantitatively monitor product attributes
- Further work needed to evaluate potential reduction in release testing panels

Potential benefit of MAM



Potential Attributes for MAM
Oxidation
Deamidation
Glycation
Glycosylation
C-terminal Lys Clipping
Sequence Variants

siRNA Structure and Chemistry Method on Specification





Method	Attribute
Non-denaturing-Ion Pair Reverse Phase by LC-UV	Duplex Purity
Denaturing-Ion Pair Reverse Phase by LC-UV	Single Strand Purity, GalNAc impurities
Denaturing-Anion Exchange by LC-UV	Single Strand Purity, PS to PO conversion
Non-denaturing-Ion Pair Reverse Phase by LC-MS	Identity

Denaturing AEX by LC-UV

Single Strand Purity



Integrated UV chromatogram




Challenges:

- Integrate entire UV chromatogram
- Report all impurities > 0.05% area
- Monitor new impurity not present in control

Alternative Format ID by LC-MS


siRNA Identification



Oligo MAM by ND-IPRP-UV-MS

Mass spec compatible method

Channel name: TUV 260



Quantification-%Purity

	% Purity		
	Control	Oxidized	Heat
Total Purity	95.46%	ND	52.64%
Loss of GalNacs	1.01% (Impurity loss of GalNac: 0.9-1.2%)	8.86%	37.86%
PS to PO	1.63%	80.73%	8.95%
Loss of nucleotides	1.32%	7.17%	ND


ND = Not Detected

% Purity=(Response/Total Response)*100

Total Purity=% Purity (Sense)+% Purity (Antisense)

Summary

- Established workflow for routine oligonucleotide ID test method
- Proof of Concept LCMS oligo MAM method to identify and quantify known impurities
- Further optimization on coeluting impurities and characterization of unknown impurities by MS



Acknowledgement

QCAS MS Development Team

Rachel Mullen

Stacey Helming

Seamus O'Connor

Paul Bigwarfe

Waters Team

Questions?

