




Tackling the Polysorbate Degradation Challenge: Control Strategies and Case Studies


Andreas Zerr
Head of Operations

WCBP 2026
29 Jan 2026

Polysorbates – complex mixtures of compounds



PS80



PS20


$$w+x+y+z=20$$

Fatty acid ester	PS20	PS80
Caproic $\text{CH}_3(\text{CH}_2)_4\text{COOH}$	$\leq 1\%$	-
Caprylic $\text{CH}_3(\text{CH}_2)_6\text{COOH}$	$\leq 10\%$	-
Capric $\text{CH}_3(\text{CH}_2)_8\text{COOH}$	$\leq 10\%$	-
Lauric $\text{CH}_3(\text{CH}_2)_{10}\text{COOH}$	40-60%	-
Myristic $\text{CH}_3(\text{CH}_2)_{12}\text{COOH}$	14-25%	$\leq 5\%$
Palmitic $\text{CH}_3(\text{CH}_2)_{14}\text{COOH}$	7-15%	$\leq 16\%$
Palmitoleic $\text{CH}_3(\text{CH}_2)_5\text{CH}=\text{CH}(\text{CH}_2)_7\text{COOH}$	-	$\leq 8\%$
Stearic $\text{CH}_3(\text{CH}_2)_{16}\text{COOH}$	$\leq 7\%$	$\leq 6\%$
Oleic $\text{CH}_3(\text{CH}_2)_7\text{CH}=\text{CH}(\text{CH}_2)_7\text{COOH}$	$\leq 11\%$	58-85%
Linoleic $\text{CH}_3(\text{CH}_2)_3(\text{CH}_2\text{CH}=\text{CH})_2(\text{CH}_2)_7\text{COOH}$	$\leq 3\%$	$\leq 18\%$
Linolenic $\text{CH}_3(\text{CH}_2\text{CH}=\text{CH})_3(\text{CH}_2)_7\text{COOH}$	-	$\leq 4\%$ <i>European Pharmacopoeia</i>

- PS are complex and heterogeneous mixtures (synthesis uses precursors from natural products)
- PS manufacturing processes can vary


Polysorbates are critical excipients for biologics

- Polysorbates are critical excipients, widely used in biologics - highly effective in protecting proteins against interfacial stress
- 56% of all marketed biotherapeutics and ~93% of all marketed mAbs contain a surfactant



Khan TA, et al, E J Pharm Biopharm, (97A), 2015, 60-67

- Polysorbates can degrade and present liabilities too:
 - Liberation of FFA → Particles → **Incompliant drug product**
 - Loss of protection against interfacial stress → **Unstable drug product / reduced shelf life**



● OXIDATION

- Presence of peroxides in the PS raw material
- Exposure to light
- Exposure to oxygen (air)
- Exposure to transition metals


● Mitigation Strategies:

- Raw material qualification, storage and handling best practices

- Ascertaining the degradation mechanism is of critical importance to implement the correct mitigation measures

● HYDROLYSIS (ENZYMATIC)

- Traces of contaminating lipases from the host cell (CAUTION: common occurrence!)



- Downstream process improvements to eliminate co-purifying enzymes

1. Reduction of PS Content


- Reduction in protection against interfacial stress
 - Instability on storage / transportation
 - Instability during processing

2. Oxidative damage of the API

(in the case of oxidative degradation)

3. Liberation of poorly soluble free fatty acids (FFA)

- Particle formation
- Potential co-precipitation of protein in biopharmaceuticals



Potential impact on

SAFETY
QUALITY
STABILITY
SHELF LIFE
COMPLIANCE



Analytical Methods for Polysorbate Characterization and Quantification

Analytical Toolbox – PS Quantification

- Quantification of Polysorbate is typically achieved by one of the following (single peak) methods:
 - Fluorescent Micelle Assay (FMA)
 - HPLC(MAX)-ELSD/ CAD
- Quantification methods need to be stability-indicating
- Quantification assays have different selectivities:
 - FMA more sensitive to oxidative degradation
 - HPLC-ELSD/ CAD - more sensitive to lipolytic degradation




Figure 6.4 Polysorbate stability indication with HPLC-FMA versus HPLC-ELSD: Analysis of 0.04% (w/v) PS80 subjected to hydrolytic stress at different lipase concentrations (top) and oxidative stress at elevated temperatures for different durations (bottom).




Figure 6.3 HPLC-ELSD applying a step gradient: Chromatogram of a 0.02% (w/v) PS80 solution.







Figure 6.2 HPLC-FMA: Chromatogram of a 0.02% (w/v) PS80 solution.

Jahn, Ch 6, Ed(s): Koulov AV, Singh SK, *Surfactants in Biopharmaceutical Development*, Academic Press, 2023, 119-136


- Polysorbate profiling by reversed phase shallow gradient HPLC (ELSD/ CAD/ MS - detection)




- Proteomic approaches to ID HCP lipases/ esterases
- Enzyme activity assays
 - PS degradation mechanism diagnosis
 - Downstream process development support




- Visible particle isolation and characterization by FTIR/ Raman microspectroscopy



- FFA target quantification by LC-MS



- Discrimination of FFA Particles by AI-aided Flow Imaging Microscopy



- PS quantification
- Polysorbate Profiling
- Lipase activity assay
- Particle isolation and ID



THE “MUST-HAVE” ANALYTICAL TOOLBOX

to ensure adequate characterization
(product quality and process control)

CONTROL STRATEGY COMPONENT	ANALYTICAL TESTS NEEDED
Raw materials testing/ qualification	<ul style="list-style-type: none">● Compendial tests● Additional characterization (profiling, FFAs)
Manufacturing process development (DS and DP), validation and in-process control	<ul style="list-style-type: none">● Quantitation (PS content)● Additional characterization (profiling, lipase assays) - DS, DP
Product characterization throughout development: Stability testing - PS quality and quantity Particles	<ul style="list-style-type: none">● Quantitation (PS content)● Additional characterization (profiling, FFAs, lipase assays, particle characterization)
Control system: Specifications (release, stability)	<ul style="list-style-type: none">● Quantitation (PS content; stability-indicating method)

clearsolutions
LABORATORIES

CASE STUDIES



CASE STUDY 1

Detection of Particles

CASE STUDY I: Appearance of Free Fatty Acid (FFA) Particles

Observation:

Visible particles
on stability



Question:


Particle composition → FFA only or
FFA + API?

Solution:

1. Particle
isolation



2. Particle
ID



- **Think ahead** (*better safe than sorry*):
 - Always **monitor PS levels on stability** (every clinical batch).
Note: do you have the right method?
 - **SvP increase** typically precedes the appearance of Visible Particles. Measuring SvPs with orthogonal methods (FIM) is a good practice
 - Consider using an **in-line administration filter** whenever possible in early clinical development
Note: Not a solution, only a proactive risk mitigation measure
- Ensure **detailed particle characterization** (stability time-course) to exclude potential protein aggregation
- Consider timing for clinical safety assessment
- Particle isolation & characterization can be technically complex - develop internal procedures OR consider a competent partner




CASE STUDY 2

Quantification Method

Observation:

- No PS quantification method available OR
- Method suitability unclear



Question:

Can a robust and stability-indicating method be developed/implemented?

Solution:

1. Review of method performance
Validation/Qualification report

2. CSL Platform PS method evaluation with the product:
a) profiling, b) quantification

3. Bespoke method development

4. Method Qualification



5. Method Transfer & Troubleshooting support

- Applied Quantification method must be stability indicating
- Especially for high concentration protein formulation the API often interferes with the PS quantification method
- Protein removal e.g. by protein precipitation could be necessary
- A robust qualification/ validation data set needs to demonstrate that the method is fit for purpose




CASE STUDY 3

Reduction of PS Content on Stability

CASE STUDY III: Elucidation of PS Degradation Mechanism

OBSERVATION:

Significant reduction of PS content on stability




QUESTION:


What is the degradation mechanism?

SOLUTION:


1. PS profiling



2. Lipase activity assay



3. FFA quantification



- Understand **the extent of the issue** (relevance of analytical results; interference)
 - Always use of **well-characterized & qualified analytical methods**
 - Do you have sufficient shelf life?
- Always use **orthogonal methods** to ensure complete picture
- Always **think about the next step** - do you have the **right analytical toolbox**?
- Consider different mitigation pathways:
 - Re-formulation - amount and type of surfactant (*be aware of pitfalls!*)
 - Shorten shelf life
 - DSP optimization




CASE STUDY 4

Lipolytic PS Degradation

CASE STUDY IV: Mitigation of Lipolytic PS Degradation

OBSERVATION:

- PS rapidly degraded
- Presence of lipolytic activity






QUESTION:

Can DSP be optimized to remove HC lipases?

SOLUTION:

Protein A optimization




1. Enzyme activity assay IPC testing to support DSP optimization

2. Select low-enzymatic-activity process version

- Consider optimizing different DSP steps
 - Protein A, HIC, MM typically give the highest probability of success
- Often multiple enzymes are present in DS (trace amounts)
- Balance removal of lipases with other attributes
- It does not have to be perfect
- Apply all characterization tools to stability assessments of the new process version - PS profiling, PS content, SvP and VP characterization

- Polysorbate degradation can present **significant challenges** to biologics development and commercial programs, with **potentially severe impact on safety and quality**
- Setting up a **comprehensive control strategy** requires:
 - An innovative, fit-for purpose analytical toolbox
 - **Use the right tool for the right job** → each question requires dedicated analytical tools

Wuchner, K et al., Do we worry too much about polysorbate degradation? An industry-wide perspective with real-life case studies, J Pharm Sci, 2026 (in press)



 expertise. insight. agility.

Thank you!