



A JSR Life Sciences Company



# Regulatory-Ready Characterization: From IND to BLA

Kelly Donovan  
Director, Characterization Services  
KBI Biopharma - Louisville, CO



# Global CDMO Providing Innovative Biologic Solutions

KBI Biopharma drives breakthroughs in biopharmaceutical development and manufacturing that positively impact the lives of patients worldwide.

**30+**  
years of  
experience

**1400+**  
client partners

**300+**  
drug candidates



**Mammalian**  
Delivering Innovative  
CLD to Vial  
Solutions



**Analytical**  
Solving Analytical  
Challenges at  
Every Stage



**Microbial**  
Elevating Microbial  
Expression and  
Manufacturing  
Efficiencies



# Product Characterization Overview

- Development Stage Expectations:
  - Investigational New Drug (IND)
  - Investigational New Drug Amendment (INDa)
  - Biologics License Application (BLA)
- ADC case study
- Best practices and potential risks



# Why Product Characterization Matters

3.2.S.3.1 Elucidation of Structure | Manufacturing Controls (CMC) | Common Technical Document (CTD)

Ultimate goal: Patient safety

## Data Collection



Primary & Higher Order Structure




Biological Activity



Variant Characterization



Degradation Pathways



## Enables

**QTPP**

Quality Target Product Profile

**CQAs**

Critical Quality Attributes

**Control**

Product and Process


**Quality**

Quality by Design

# Biotherapeutic Development Phases



# Biotherapeutic Development Phases



# From IND to BLA: Increasing Analytical Expectations

## IND EoS

- **Purpose:** Initial Characterization Assessment for First in Human (FIH) trials
- **Material:** Tox/Clinical Lot
- **Methods:** Platform
- **Duration:** 6 weeks

## INDa EoS

- **Purpose:** Deeper Product Understanding in Phase II/III
- **Material:** Phase I/II DS and DP
- **Methods:** Optimized for product
- **Duration:** 3 – 5 months

## BLA EoS

- **Purpose:** Full Package for Commercial Approval
- **Material:** Commercial Scale/PPQ Lot
- **Methods:** Qualified, product specific
- **Duration:** 9 - 12 months

Increasing Complexity

# IND: Early-Stage Structural Characterization


- Fast characterization focused on patient safety
- Approval required for Phase 1 (FIH) trials

Category	Method	Purpose
Primary structure	Intact mass	Mass measurements
	Reduced peptide mapping	Post Translational Modifications (PTMs)
Higher order structure	FTIR/Far UV CD spectroscopy	Secondary structure (i.e. alpha helices, beta sheets)
	Near UV CD spectroscopy	Tertiary structure
Biological Activity	Binding or Potency method	Structure function

Methods presented are based on mAb therapeutics. More complex biologics may utilize additional characterization techniques.

# INDa: When Characterization Becomes Critical

## INDa: Required for significant process changes



### Common Triggers:

- Scale up (i.e., lab → pilot → GMP)
- New manufacturing site or CDMO
- Raw material changes (e.g., cell culture media)

### • Requirements:

- Comparability data using product-specific methods
- Multiple lots to demonstrate no new safety risks

# INDa: When Characterization Becomes Critical



## Primary Structure

- Intact mass
- Peptide mapping
- Disulfide mapping
- Glycan profiling



## Size Variants

- SEC
- SEC-MALS
- CE-SDS



## Biological Activity

- FcRn binding
- Potency assay



## Higher-Order Structure

- FTIR/Far UV CD spectroscopy
- Near UV spectroscopy
- DSC



## Charge Variants


- IEX
- cIEF



## Particle Assessment

- MFI (Drug Product)

# BLA: Complete Characterization Package




# Characterization Requirements: Side-by-Side View

Increasing Complexity 


Feature	IND	INDa	BLA
<b>Primary Structure</b>	<ul style="list-style-type: none"> <li>Intact mass</li> <li>Reduced peptide map</li> </ul>	<ul style="list-style-type: none"> <li>Intact mass</li> <li>Reduced peptide map</li> <li>Disulfide map</li> <li>Glycan map(s)</li> </ul>	<ul style="list-style-type: none"> <li>Intact mass</li> <li>Reduced peptide map(s)</li> <li>Disulfide map</li> <li>Glycan map(s)</li> </ul>
<b>High Order Structure</b>	<ul style="list-style-type: none"> <li>FTIR/Far UV CD</li> <li>Near UV CD</li> </ul>	<ul style="list-style-type: none"> <li>FTIR/Far UV CD</li> <li>Near UV CD</li> <li>DSC</li> </ul>	<ul style="list-style-type: none"> <li>FTIR/Far UV CD</li> <li>Near UV CD</li> <li>DSC</li> </ul>
<b>Impurity Characterization</b>	Not required	Charge characterization of peak groups	Charge & size characterization of individual peaks
<b>Biological Activity</b>	Binding or potency assay	Binding & potency assays	Binding & potency assays
<b>Degradation Pathways</b>	Limited for early understanding of molecule stability	Limited to support other studies: comparability, formulation dev, method qualifications	Full degradation of key stress conditions: heat, pH, oxidation, light, agitation
<b>Method Qualification</b>	Not required	Limited, focus on method optimization for product	Qualification for characterization methods

# Project Timelines: How to be Efficient




## Couple RS Characterization and/or Comparability with Product Characterization

# Project Timelines: How to be Efficient



## Couple RS Characterization and/or Comparability with Product Characterization

# Project Timelines: How to be Efficient



## Couple RS Characterization and/or Comparability with Product Characterization



A JSR Life Sciences Company


# Case Study: ADC Characterization




Analytical

# ADC Structure Includes Multiple Components for Characterization

CQAs for **antibody**, **drug-linker**, & combined conjugated product




# Commonly Known ADC CQAs




- Drug-to-Antibody Ration (DAR)

# Common ADC CQAs




- Drug-to-Antibody Ration (DAR)
- Conjugation Site

# Common ADC CQAs




- Drug-to-Antibody Ration (DAR)
- Conjugation Site
- Free Drug

# Common ADC CQAs




- Drug-to-Antibody Ration (DAR)
- Conjugation Site
- Free Drug
- Linker Stability

# Common ADC CQAs



- Drug-to-Antibody Ration (DAR)
- Conjugation Site
- Free Drug
- Linker Stability
- Aggregation/Size Variants

# Common ADC CQAs



- Drug-to-Antibody Ration (DAR)
- Conjugation Site
- Free Drug
- Linker Stability
- Aggregation/Size Variants
- Potency

# IND: Fast Assessment for FIH Safety



## DAR

- HIC



## Free Drug

- Intact mass



## Size Variants

- SEC
- CE-SDS



## Conjugation Site

- Peptide mapping



## Linker Stability

- HIC – changes in DAR levels



## Potency

- Binding assay
- Cytotoxicity assay

Testing completed on using platform methods

# INDa: Demonstrate Comparability For Site Change



## DAR

- HIC
- Intact mass



## Free Drug

- Intact mass
- RP-HPLC



## Size Variants

- SEC
- CE-SDS
- SEC-MALS



## Conjugation Site

- Peptide mapping



## Linker Stability

- HIC – changes in DAR levels



## Potency

- Binding assay
- Cytotoxicity assay
- Cell based bioassay

**Testing completed on using product specific methods**

# BLA: Complete Characterization Package for Commercial Approval



## DAR

- HIC
- Intact mass
- Fractionation of hydrophobic variants



## Free Drug

- Intact mass
- RP-HPLC



## Size Variants

- SEC
- CE-SDS
- SEC-MALS
- Fractionation of size variants
- SV-AUC



## Conjugation Site

- Peptide mapping
- Intact mass



## Linker Stability

- HIC – changes in DAR levels
- Degradation studies



## Potency

- Binding assay
- Cytotoxicity assay
- Cell based bioassay

Testing completed on qualified characterization methods

# Avoidable Mistakes That Delay Filings



## **Inadequate Qualification of Characterization Methods**

- Not qualifying LC-MS and HOS methods
- Lack of method performance understanding at late-stage

## **Incomplete Characterization Studies**

- Assessing only size OR charge variants (not both)
- Skipping hydrophobic variants when detected

## **Delayed Characterization Timelines**

- Starting characterization too late in development
- Discovering critical issues at BLA stage

## **Insufficient Comparability Data**

- Incomplete methods to demonstrate impact on CQAs
- Not enough lots for stage of comparability

# Regulatory Expectations - Where the Bar Is Heading

## Today

Basic IND Data

Late CQA  
Identification

Manual Analysis

Raising  
the Bar

## Future

Robust IND Data  
(Complex Modalities)

Early CQA Identification  
Risk-based Strategy

AI Enabled Modeling

# Bringing It All Together

- Product characterization is progressive – start light, end deep
- IND: Fast with a focus on patient safety
- INDa: Comparability studies with greater product understanding
- BLA: Complete package
- Regulatory success = scientific rigor + staged planning



# Virtual Poster Presentation

65

## LC-MS Characterization of an Unexpected C-Terminal HC Extension in a Therapeutic mAb

Topic: WCBP Virtual Poster

### LC-MS characterization of an unexpected C-terminal HC extension in a Therapeutic mAb

Ben Zigterman



#### Introduction

To ensure safety and efficacy therapeutic monoclonal antibodies (mAbs) require extensive molecular characterization, including identification and characterization of size variants. SDS-PAGE of an IgG1 drug candidate revealed an unexpected HMW band under reducing conditions.

Orthogonal IEF and SEC did not observe this anomaly, suggesting it comprised a HC size variant. Accordingly, LC-MS was conducted to distinguish among likely modifications (e.g. contaminants, glycovariants, and frameshift mutants).

#### Results

**Mass Analysis Reveals Species +2448 Da and +2520 Da Greater than Expected, with Shift Localized to Fc**

Species	Major Form	Anomaly	A
IgG (GFP-FC)	54399.4	543524.1	2518
HC	48847.5	51187.8	2520
FC (GFP)	25204.8	27724.6	2520

*Note: Data for light chain and IgG separated FC were identical between the 2 lots and are not shown here.*

*Unconvoluted masses of intact IgG (A), deglycosylated heavy chain (B), and IgG1C digested Fc (C) shown above for each lot, showing the presence of mass variant exclusively in lot 2 and localizing the variant to the Fc.*

**The variant Fc resolved chromatographically, and the spectra of the variant peak is shown in panel (D).**

#### Results Continued

**Denovo Sequencing of C-terminal Extensions**

*Note: Leucine and Isoleucine are indistinguishable by this methodology so any 'L' above could be explained alternatively as an 'I'. MS<sup>n</sup> data were sequenced de novo then compared the theoretical sequence from the expression vector, revealing that unknown species comprised C terminal extensions caused by alternative splicing of the vector sequence.*

**Alternative Splicing Yields C-Terminal Extensions**



A JSR Life Sciences Company

# Science Driven. Customer Focused.

A Global CDMO Providing Innovative  
Biologic Solutions



Kelly Donovan



[kdonovan@kbibiopharma.com](mailto:kdonovan@kbibiopharma.com)



+1 720-845-6002



[kbibiopharma.com](http://kbibiopharma.com)