

Quantitative Cell Characterization Using Laser Force Cytology (LFC): The Future of Precision Analytics for Improved Manufacturing Success

Colin Hebert, PhD Senior Vice President, Scientific and Business Operations

CASSS WCBP 2025 January 28, 2025

#### LASER FORCE CYTOLOGY™

## **Biomanufacturing Advanced Therapies is Complicated**

Starting materials are complex, variable, and difficult to characterize, reducing process and product consistency. Delayed and imprecise analytics result in lost process and production efficiency, delaying delivery of life saving treatments. Large biological variability throughout manufacturing processes results in batch-to-batch variability and potential for OOS/OOT and batch failure events.

Current Analytics

Slow 🕂

•

•

Robust cellular characterization throughout a process paves the way for predictive insights, rapid optimization and adaptive manufacturing

LFC Technology Background Vaccines and Viral Vectors

Stem Cell Therapies CAR T Cell Therapies



2 LumaCyte Laser Force Cytology © LumaCyte 2025 Proprietary & Confidential



## e Real-time Label-free Single Cell Analysis Provides Solutions

- Label-free cellular PAT analytics reduce the need for antibodies, provide unbiased measurements and enable the sensitive measurement of cellular changes
- Accurate, precise and non-subjective multivariate data comprehensively characterize complex starting materials and products
- **Real-time analytics increase process knowledge**, enable improved process controls and optimization, and maximize product quality, yield, and shelf life
- **Quantitative precision analytics allows biological evidence** to be carried forward from early process development, production and QC
- Reduced labor and waste significantly lowering costs and accelerates time to patient







Label-free, single cell analysis based on intrinsic biochemical & biophysical properties using a balance of optical and hydrodynamic forces in a microfluidic channel



Virally infected Vero cells

Laser (Optical Force) → Velocity ∝ F<sub>Optical</sub> <------ Fluid Flow (Drag Force)

- An optical force is generated when a laser beam reflects and refracts through a cell
- Laser Force Cytology<sup>™</sup> measures velocity (optical force) and other parameters to detect subtle phenotypic changes in cells, rapidly measuring quantitative early indicators of cellular response to viral infection, activation, transfection and differentiation
- Applications in **cell and gene therapy, vaccine development, and biomanufacturing**



#### HOW DOES LASER FORCE CYTOLOGY<sup>™</sup> CONNECT TO BIOLOGY?





Potency

**Cellular Fingerprinting** 



Reilly McCracken<sup>1</sup>, Noor Al-Nazal<sup>1</sup>, Travis Whitmer<sup>1</sup>, Sijia Yi<sup>1</sup>, James M. Wagner<sup>2</sup>, Colin G. Hebert<sup>3</sup>, Matthew J. Lowry <sup>3</sup>, Peter R. Hayes <sup>4</sup><sup>(0)</sup>, James W. Schneider <sup>4</sup>, Todd M. Przybycien <sup>5</sup> and Malini Mukherjee <sup>6,\*</sup> Global Vaccines and Biologics Commercialization, Merck & Co., Inc., West Point, PA 19486, USA Vaccines Process Development, Merck & Co., Inc., West Point, PA 19486, USA LumaCyte, Inc., 1145 River Road, Suite 16, Charlottesville, VA 22901, USA Department of Chemical Engineering, Carnegie Mellon University, Pittsburgh, PA 15213, USA Department of Chemical and Biological Engineering, Rensselaer Polytechnic Institute, Troy, NY 12180, USA Analytical Research & Development, Merck & Co., Inc., West Point, PA 19486, USA Correspondence: malini.mukheriee@merck.com to prevent infectious diseases are given to target the body's innate and Volume 36, Issue 26, 18 June 2018 ISSN 0264-410X ms. In most cases, the potency of a live virus vaccine (LVV) is the most of efficacy, though in some cases the quantity of surface antigen on the virus uality attribute. Existing methods to measure the potency of viruses include ays, both of which have very long lead times and cannot provide real time ality of the vaccine during large-scale manufacturing. Here, we report the yte's Radiance Laser Force Cytology platform as a new way to measure the accine am biomanufacturing process in real time and compare this to traditional o assess this new platform as a way to detect adventitious agents, which tion for the release of commercial vaccines. In both applications, we report pedited and relevant potency information with strong correlation to release ether, our data propose the application of Laser Force Cytology as a valuable Science\_ 12 DECEMBER 2018 Translational Medicine

Transfection

**Donor variability** 

Transduction

Cell health/viability

Electroporation

Viral Infection and

Neutralization

#### Cryopreservation **Formulation/Processing**

Hebert et al. 2011. Anal. Chem. 83, 5666-5672

MAAAS

MDPI

### **Cadiance** Broad Applicability Delivering Robust Results *Cell, Virus and GOI Independent*

Systems Tested with Laser Force Cytology to Date

#### Mammalian Cells and Viruses

(partial list due to client confidentiality):

- Vero, MDCK, HEK-293, MRC-5, MARC 145, CRFK, A549, HL60, U937, SK-MEL-28, HCC2429, HEp-2, CHO, Primary Human Blood components including PBMCs, T-cells and RBCs, HeLa, NALM-6, iPSCs, hMSCs, cardiomyocytes
- Influenza, VSV, Zika Virus, VZV, PRRSV, Swine Influenza Virus, Pseudorabies Virus, Feline Calicivirus, Adenovirus, HSV-1 (T-Vec), RSV, Mouse Minute Virus, retrovirus, lentivirus, AAV, human coronavirus, Measles virus

#### Insect Cells

 Baculovirus infection of Sf9 cells for vaccine, VLP, protein, and AAV production



Reproducible measurements of the same samples over a 50+ day period demonstrate instrument consistency



Average Velocity: 2393 ± 11 Uninfected Vero cells (0.46% CV) 2244 ± 14 VSV Infected Vero cells (0.36% CV)



## **Regulatory Support for Novel PAT Analytics**





Peter Marks, MD, PhD

**Director of the Center** 

for Biologics Evaluation

and Research, FDA

#### **GEN** Webinars

**Rapid Bioprocess Analytics Needed to Drive Improvements** in Product Consistency and Quality

**Fixing The Potholes in** 





Approach for Establishing Novel and Fit-for-Purpose Cell Viability Methods to Support Cell Manufacturing Process Monitoring MATERIAL MEASUREMENT LABORATORY

armona Sarkart, Charles H. Camp Jr.\*, Matthew Lowry\*, Laura Pierce\*, David Varisco\*, Colin Hebert\*,, Zachary Evans\*, Sean Hart\*

An optical force is generated when a laser



versells (e.g. 89 cells) can lead to high

result in classifiers with poor accuracy



Test Samples with Systematically Varied Cell Health and Cell Growth Profiles rist cells were cultured in different media conditions, resulting in a series of test samples, with systematically varying cell health conditions (visible cell state





#### **Analytical Procedures for Viral Vectored Vaccine Quality**

Draft guidelines

Accelerating product development through a common understanding of quality





## Vaccine and Viral Vectors Cell Health, Potency, and Cell-Based PAT

Radiance'



### Cell Viable State Monitoring Predictive Insights for Consistent Process Performance



- Radiance<sup>®</sup> can rapidly monitor cell health conditions and characterize changes in cell populations to ensure consistent process performance
- This study compares Vero cells frozen with different methods to the standard protocol, cells were then analyzed on Radiance<sup>®</sup> immediately post thaw
- Significant changes are seen in Radiance<sup>®</sup> Optical Force Index for samples without DMSO (Sample 3) and rapid (uncontrolled) freezing (Sample 4)

Sample	Condition		
1	Standard/Same as reference		
2	Higher FBS concentration		
3	No DMSO		
4	Uncontrolled freezing		



#### Measuring H<sub>2</sub>O<sub>2</sub> Induced Oxidative Stress Improving Cell Health Measurement with LFC



- Vero cells treated with various levels of H<sub>2</sub>O<sub>2</sub> to simulate oxidative stress felt during bioproduction
- Radiance detects changes in velocity and eccentricity in all conditions compared to control, with no statistical difference seen in viability between untreated cells and up to 300 µM H<sub>2</sub>O<sub>2</sub>
- Caspase activity assay demonstrates that this LFC metric corresponds more closely with apoptosis, providing additional information beyond trypan blue viability



#### Reduced Viral Production Due to Oxidative Stress Fit for Purpose Cell Health Measurements



24h Virus		
Production	Sample	VSV Production (TCID50/mL)
	Control (0 µM H <sub>2</sub> O <sub>2)</sub>	1.31E10
	450 µM H <sub>2</sub> O <sub>2</sub>	1.25E10
	750 μM H <sub>2</sub> O <sub>2</sub>	4.22E9
	1500 µM H <sub>2</sub> O <sub>2</sub>	1.31E9

- Vero cells treated with various levels of H<sub>2</sub>O<sub>2</sub> to simulate oxidative stress felt during bioproduction
- Stressed and untreated cells frozen, thawed, and infected with vesicular stomatitis virus (VSV) ~6h post thaw
- Effects of the stress are visible at each stage of the process, including as a greater than 50% reduction in viral production at the 750 μM condition compared to unstressed cells





### **Compressing Time to Results** *Rapid Viral Analysis with Radiance® vs TCID50*

- Faster infectivity / titer results:
  1 day versus 10 days (adenovirus)
- Correlates with plaque assay or TCID50 results
- Absolute infectivity measurements without correlation
- Direct results from bioreactor with no further incubation (answers in 5 minutes)
- Radiance<sup>®</sup> delivers higher precision, removing subjectivity to streamline assay transfer and comparability



## **Rapid Viral Analysis** LumaCyte Vesicular Stomatitis Virus (VSV) Infectivity in Vero Cells



Radiance<sup>®</sup> detects subtle phenotypic changes in virally infected cells, rapidly measuring quantitative early indicators of cellular response to viral infection.



## **Calculating Viral Titer with Radiance**<sup>®</sup> *VSV in Vero Cells*



- Normal titer process with TCID50 takes 72h post infection
- Radiance takes only 16 hours for a MOI range of 0.0001 to 1
- Infection metric can be used directly or correlated to TCID50



Hebert et al. 2018. Vaccine, 36, 6061–6069







Stock TCID50/mL: 1.56E+8

- Composite calibration curve developed between Radiance<sup>®</sup> data and MOI across 5 independent experiments
- Separate unknown samples were predicted using the calibration curve with an average log10 difference of 0.06
- ✓ Demonstrates rapid analysis in <24h post infection with high accuracy and precision vs 72h for TCID50



### **Readiance** Adenovirus Absolute Titer Analysis 48 Hours Post Infection



Plaque Assay Titer: 3.61E+09

	Radiance Absolute Titer	Log10 Difference
Ехр А	3.51E+09	-0.012
Exp B	4.03E+09	0.048
Exp C	5.03E+09	0.144
Average	4.19E+09	0.065
St Dev (CV%)	7.73E+08 (18%)	

- Absolute titer algorithm can be used to calculate viral infectivity without correlation
- Results in table show absolute average log<sub>10</sub> difference of 0.065 (14%) when comparing absolute titer and plaque assay results
- Absolute titer results in 48h post infection versus 7+ days for plaque assay



### VSV Physical Count vs Infectious Titer Results Summary



Treatment	<u>NanoSight Pro:</u> Total Particle Count (particles/mL)	<u>Radiance:</u> Infectious Viral Titer (TCID50/mL)	Log <sub>10</sub> Reduction from Untreated
Untreated	2.51E9±2.80E8	3.31E9±1.03E9	N/A
Freeze Thaw	1.81E9±1.80E8	1.61E9±6.02E8	0.31
Low pH	1.74E9±4.70E8	6.83E6±4.95E6	2.69
High Temp	2.78E9±3.60E8	3.40E6	2.99

Webinar Available on Demand

- Vesicular stomatitis virus (VSV) samples subjected to a variety of insults that may affect product quality
- Effects of the insults quantified by viral titration using Laser Force Cytology (Radiance) and particle counting using Nanoparticle Tracking Analysis (NanoSight Pro)

Total Particle Count

Infectious Viral Titer

No appreciable changes in total particle count, but Low pH and High Temperature show a significant reduction in viral titer, demonstrating the importance of measuring functional viral titer





- Under reduced serum conditions, Gibco Production Medium No. 1 and No. 2 show improved viral titer compared to Advanced MEM
- Highest viral titer achieved at 24h time point (~14X higher than standard conditions)
- Rapid infectivity assay enables screening of wide process design space



## **Radiance** AAV Multivariate Transduction Curve



• Strong correlation developed between Radiance data and TCID50/mL titer

 Accurate and precise unknown prediction within 0.045log<sub>10</sub> (9.4%) of the known stock titer with a CV of 14% across 4 unknows





## **Radiance** Process Analytics and Production

Rapid and Precise Process Analytics Enable:

- Real-time process optimization
- Increased process knowledge and understanding of process variability
- Real-time contamination monitoring and early mitigation of potential batch failures
- Significant increases in product yields





- Example viral process optimization experiment to determine the optimal multiplicity of infection (MOI)
- Relative response of Radiance infection metric can be monitored in real-time to rapidly screen process conditions without having to perform lengthy and resource intensive assays such as plaque/TCID50, PCR, or ELISA



## Live Virus Vaccine Production Monitoring Estimated Potency Correlation



Experiment	Time Point (h)	Radiance Infection Metric (% OFI>55 s <sup>-1</sup> )	Estimated Potency (TCID50/viable cell)	Calculated Potency (TCID50/viable cell)	Log <sub>10</sub> Difference
	72	11.53	1.36	1.55	0.056
Exp 1	96	23.52	2.97	2.54	-0.068
	120	40.92	Estimated Potency (TCID50/viable cell)      Calculated Potency (TCID50/viable cell)        1.36      1.55        2.97      2.54        6.81      5.22        1.38      1.42        1.98      2.45        3.79      5.59        10.11      8.97        3.06      3.08        4.85      3.28        5.04      5.34	-0.116	
	72	9.49	1.38	1.42	0.015
Evp 2	96	22.60	1.98	2.45	0.093
Exp 2	120	42.57	3.79	5.59	0.169
	144	54.02	10.11	8.97	-0.052
	72	28.12	3.06	3.08	0.002
Even 2	96	29.69	4.85	3.28	-0.169
Exb 2	120	41.48	5.04	5.34	0.025
	144	65.92	13.22	14.67	0.045
				Absolute Average	0.074

McCracken et al. 2022. Vaccines 10(10), 1589

- O Work conducted by Merck to monitor real-time viral infectivity from a bioreactor
- Strong correlation demonstrated between estimated potency measurements and Radiance data
- Calibration curve can be used to accurately calculate the estimated potency throughput the production process and provide real time feedback to improve the speed of process development and manufacturing consistency
- Additional data in publication demonstrates monitoring of another virus, real-time contamination monitoring, and rapid analytics (TCID50 replacement) for measles virus





### Stem Cell Applications Label free characterization of differentiation

Radiance'

AAV



### **Cadiance** Label-free iPSC Differentiation Monitoring LumaCyte LFC Metric Changes



- Induced pluripotent stem cell (iPSC) samples were directed into three different lineages (endoderm, mesoderm, and ectoderm) and compared to undifferentiated controls as well as non-directed/spontaneous samples
- Multivariate analysis can be used to visualize changes seen in multiple parameters



### **Readiance** Label-free iPSC Differentiation Monitoring *Principal Component Analysis*



PCA: Multivariate dimensionality reduction technique to easily visualize multidimensional data

- Principal component analysis (PCA) uses multiple Radiance data parameters to assess cell populations each point represents a population of cells
- Radiance can be used to monitor iPSC and other cell differentiation pathways in a label-free manner, allowing for unbiased characterization of cells without the limitations of antibodies that only indicate the presence of a specific marker or protein
- Allows for quantitative and universal characterization of cell differentiation





## **Radiance** Label-free hMSCs Differentiation Monitoring



 LFC single cell data can be used to monitor hMSC differentiation by tracking population shifts and identifying specific sub-populations of interest

# Readiance<sup>®</sup> Label-free hMSCs Differentiation Monitoring



Study to monitor the expansion and Ο differentiation of human bone marrow derived mesenchymal stem cells into osteoblasts, adipocytes and chondrocytes

Radiance can be used to monitor MSCs  $\checkmark$ differentiation in a label-free manner at different time points enabling a rapid and label-free indication of commitment towards a specific lineage











27 LumaCyte Laser Force Cytology © LumaCyte 2025 Proprietary & Confidential

#### **Cadiance** LumaCyte Monitoring hMSC Changes due to Passage Experimental Summary

- Study conducted to investigate how hMSCs differentiation changes as passage number increases
- Undifferentiated cell samples analyzed using Radiance at passage numbers 5, 8, and 16
- At each of the three passages, cells were induced to differentiate into adipocytes and stained 14 days post differentiation
- <u>Objective</u>: determine if differences in undifferentiated cells can be measured using LFC that could be used to predict success or failure of differentiation



### **Cadiance** Monitoring hMSC Changes due to Passage LumaCyte LFC Population Data



- Changes seen in multiple LFC metrics, including both average velocity and average complexity factor as passage number increases
- Velocity is proportional to the cell's optical force, while complexity factor is a shape metric related to the cell's circularity
- As shown on the next slide, multivariate methods can be used to measure population wide changes in multiple metrics



#### **Readiance** LumaCyte Monitoring hMSC Changes due to Passage Principal Component Analysis



- O PCA plot of undifferentiated cell populations at different passage numbers shows statistical separation between samples
- LFC data could be used to predict future differentiation quality





Passage 16



#### CAR T Cell Therapy Co-culture killing, in-process phenotype, and predictive donor characterization

Radiance'



#### LFC Analytics in Manufacturing – CAR T Cell Therapy Example



O Indicates where Radiance can be used across complex cell therapy manufacturing processes

\* CAR T manufacturing process can vary between developers (4-12 days)



#### LFC Analytics in CAR T Cell Therapy Manufacturing



() Indicates where Radiance can be used across complex cell therapy manufacturing processes

\* CAR T manufacturing process can vary between developers (4-12 days)



### Process Tracking using PCA T-cell Activation



- T-cell data across 5 experiments and 3 different donors used to create process "road map" using LFC data
- Deviation in process detected outside boundary as early as 1 day post activation and continues down alternate path

Process average and
standard deviation:
$3.2 \sigma$ in PC1
3.0 σ in PC2

 $\bigcirc$ 

**Deviation: Doxorubicin treated T Cells** 



34 LumaCyte Laser Force Cytology © LumaCyte 2025 Proprietary & Confidential

## **Readiance** Process Prediction using PCA *T-cell Activation*

- Analysis of cells at early timepoints can predict future process failure
- 24 hr PCA model shows differences in cell populations affected by treatment with a chemotherapeutic agent versus untreated cells, providing near real-time feedback of a process variation



## Poseida Therapeutics using ML models to predict cell state based on LumaCyte Radiance® data





Accurate mapping of Radiance<sup>®</sup> data to cell health and cell phenotype.

Enables real-time monitoring of cell state in process

#### Model Application

Model outputs interpret/predict the result without the need to run flow cytometry

36 LumaCyte Laser Force Cytology

 $\ensuremath{\mathbb{C}}$  2025 Poseida Therapeutics, Inc. All Rights Reserved.



#### Poseida Therapeutics utilization of LFC Process Analytical Technology (PAT) and AI in developing a control strategy



- ML models exhibit excellent predictive power on cell health and phenotype across multiple in-process timepoints and donor starting material.
- Correlations between certain product attributes or process performance with input parameter or combination of input parameter(s) (e.g. PCA principal component analysis) can be made
- These correlations can be utilized for control strategy
- Results from PAT can be applied to guide manufacturing decisions at critical process steps



© 2025 Poseida Therapeutics, Inc. All Rights Reserved.

## Detection of CAR T Mediated Cell Killing Average Velocity (x optical force)



- CD19 CAR T cells were produced from healthy donors and killing efficacy was determined using a 24h co-culture assay with NALM-6 luciferase target cells at a variety of target to effector cell ratios
- Specific detection of cell killing in CN population when compared to TN and MN populations was detected as an increase in velocity
- A decrease in killing response is seen as fewer effector cells are added to the co-culture



## Variation of Cell Killing Across Donors



- Changes in LFC metrics across populations can be used to measure cell killing across donors
- Variation in target cell (NALM6) killing as a function of number of CAR T cells across 5 healthy donors after a 24h coculture period
- Highest overall killing shown in Donors B and E, especially at the lower ratios



# What if you could predict CAR T potency based upon the starting materials?





Successful CAR T production



Outcomes-based reimbursement (OBR) for high \$ therapies<sup>1</sup>

#### Unsuccessful CAR T production



- Ineffective treatment
- Wasted effort and \$
- Delays
- Negative prognosis

<u>1. Jesper Jørgensen</u>, <u>Eve Hanna</u>, and <u>Panos Kefalas</u>, "Outcomes-based reimbursement for gene therapies in practice: the experience of recently launched CAR-T cell therapies in major European countries", <u>J Mark Access Health Policy</u>. 2020; 8(1): 1715536.



#### **Readiance** LumaCyte Correlation of CAR T Potency with Cellular Starting Material



Ellipses are  $2\sigma$  boundaries

	Percent Killing (%)				
Ratio (CAR T:NALM6)	Donor A	Donor B	Donor C	Donor D	Donor E
1:1	97.32	99.68	98.69	93.02	98.56
0.5:1	95.57	99.43	96.83	87.53	97.48
0.2:1	51.59	82.14	41.54	35.23	78.16
0.1:1	0.50	67.62	0.50	5.49	54.14

- PCA model created using Radiance data collected from donor PBMC samples
- Donors B and E are the high killing donors and the others are the lower (see table)
- Radiance may be used to estimate the killing efficiency of donor T-cells prior to manufacturing of the therapy



# What if we could use Prescriptive Analytics for Adaptive Manufacturing?



"When feasible, risks to potency caused by material variability should be mitigated by designing a manufacturing process with <u>adaptive steps</u> that compensate for variations in the material."

FDA Potency Assurance for Cellular and Gene Therapy Products – Draft Guidance for Industry, December 2023, 338 – 340







The development and manufacturing of cell-based products is complex and requires precision analytics to reduce costs and increase access



Current analytical tools are often slow, laborious, and do not always predict or indicate manufacturing or clinical success, especially when characterizing cellular material



Laser Force Cytology (LFC) uses a combination of optical and fluidic forces to sensitively and rapidly measure cellular properties and responses in a label-free manner



LFC can be used throughout the development and manufacturing process, including correlation of cellular starting material quality to predict manufacturing success



LumaCyte's Radiance<sup>®</sup> instrument enables rapid assessment of cell health, phenotype, viral infectivity and cytotoxicity to facilitate rapid optimization and adaptive manufacturing







## Thank You For Your Time & Attention

#### Acknowledgements

#### LumaCyte

Mina Elahy Zach Evans Anna Kersey Matt Lowry Jonathan Turner Sandra Walton Millicent Williams

info@lumacyte.com

Merck and Co. NIST Poseida Therapeutics

LumaCyte.com