

Live the Legacy.
Protect the Future.



Development of a Reporter Cell Line and Assay for Merck's Ebola Vaccine (rVSV- Δ G-ZEBOV-GP)

Presented by Brian K. Meyer
Merck & Co., Inc., West Point, PA, USA

Bioassays 2023
April 17-19



Overview

- Recombinant Vesicular Stomatitis Virus Ebola Vaccine (rVSV- Δ G-ZEBOV-GP)
- Background and rationale for developing a cell-based reporter assay vs. plaque potency
- Reporter cell line development: JM-1
- Assay development with JM-1 and rVSV- Δ G-ZEBOV-GP
- Conclusions
- Acknowledgements
- References

Recombinant Vesiculostomatitis Virus (rVSV)-Based Ebola Vaccine

- rVSV was generated in a cell line using reverse genetics, which included the glycoprotein gene from the Ebolavirus-Zaire, Kitwit strain
- rVSV- Δ G-ZEBOV-GP is grown in cell culture followed by purification, manufacture of drug substance, and drug product
- Infectivity by plaque titration was utilized to determine final titers of the drug substance, drug product, and stability samples
- A cell-based reporter assay was investigated as an alternative to infectivity by plaque titration for rVSV- Δ G-ZEBOV-GP





Figure 1. Generation of rVSV- Δ G-ZEBOV-GP

Background and Rationale: Cell-based Reporter Assays vs. Plaque Potency

Goal: Determine if a cell-based reporter assay could be developed that would be more efficient than the plaque potency assay for a Merck rVSV-based vaccine (Ebola, rVSV- Δ G-ZEBOV-GP)

- Plaque assay
 - The industry standard for determining potency of live virus vaccines
 - Assay infectivity via plaque titration
 - Low throughput, requires significant resources, and has high variability
 - Long turnaround times can slow down vaccine development
- Cell-based reporter assays
 - Designed to be fast and with low variability
 - Acceleration of viral vaccine development is critical in the event of an outbreak



Ref 1




Figure 2. A) Plaque assay
B) Cell-based reporter assays⁴

Reporter Cell Line Development: JM-1

Goal: Develop a reporter cell line that could be used with rVSV- Δ G-ZEBOV-GP

- Decision was made to generate a cell line that would constitutively express a reporter protein in the cell, selected NanoLuc (Promega, engineered from the deep sea shrimp *Oplophorus gracilirostris*²)
- Generated a vector construct consisting of the CMV promoter driving NanoLuc. Isolated stable transfectants using G418 selection and cloned JM-1
- Hypothesis was NanoLuc would be released from the cell following infection by rVSV- Δ G-ZEBOV-GP
- Demonstrated proof-of-concept with rVSV- Δ G-ZEBOV-GP (Ebola vaccine, ERVEBO®)³




Figure 3. A) Illustration of JM-1 and release of NanoLuc following infection with rVSV- Δ G-ZEBOV-GP, B) Cell culture plate of JM-1 infected with rVSV- Δ G-ZEBOV-GP and substrate for NanoLuc

JM-1 Cell Line: Test Proof-of-Concept

Goal: Evaluate proof-of-concept that an increase in rVSV- Δ G-ZEBOV-GP would result in a concomitant increase in NanoLuc expression in the media

- Selected a condition to obtain proof-of-concept
- A range of virus was tested in a 96-well plate
- Substrate added and luminescence measured




Figure 4. JM-1 infected with an increasing titer of rVSV- Δ G-ZEBOV-GP

Conclusion: Proof-of-Concept completed, an increasing titer of rVSV- Δ G-ZEBOV-GP resulted in an increase in NanoLuc

JM-1 Cell Line: Stability-indicating

Goal: Determine if JM-1 was stability-indicating for rVSV- Δ G-ZEBOV-GP by measuring NanoLuc activity in the media

- Samples of rVSV- Δ G-ZEBOV-GP were incubated at accelerated temperatures for 0, 3, and 7 days
 - Accelerated temperatures will result in lower potency of rVSV- Δ G-ZEBOV-GP
- Samples were assayed and media measured for NanoLuc activity





Figure 5. JM-1 is stability-indicating for rVSV- Δ G-ZEBOV-GP

Conclusion: JM-1 was stability-indicating

JM-1 Cell Line: Optimize conditions for infection

Goal: Optimize cells/well and incubation time for assay development


- Evaluated 2E4, 4E4, and 6E4 JM-1 cells/well
- Infection time of 30, 48, and 66 hours, followed by luminescence measurement



30 Hours	R ²
20k cells/well	0.5518
40k cells/well	0.6573
60k cells/well	0.7232



48 Hours	R ²
20k cells/well	0.8844
40k cells/well	0.9236
60k cells/well	0.9511



66 Hours	R ²
20k cells/well	0.5820
40k cells/well	0.9092
60k cells/well	0.8766

Figure 6. Experimental results of cell number and infection time with rVSV-ΔG-ZEBOV-GP

Conclusion: 6E4 cells/well infected for 48 hours had the highest R² under these experimental conditions

JM-1 Cell Line: Evaluate correlation to plaque potency assay

Goal: Determine if there was a correlation between the JM-1 and plaque potency assays

- Prepared a sample set of rVSV-ΔG-ZEBOV-GP, submitted to the plaque potency lab, and tested with JM-1





Figure 7. A) Experimental design, B) Concordance analysis of the NanoLuc vs. plaque potency

Conclusions: Determined there was a correlation between the NanoLuc and plaque potency assays⁹

Evaluate optimum dilution for curve generation

Goal: Determine the optimized dilution to obtain a standard curve and evaluate parallelism



B)




Figure 8. A) Experimental design for rVSV-ΔG-ZEBOV-GP positive control and reference and four parameter fit with 1:4 and 1:5 dilutions; B) Common slope and Individual fits models

Conclusion: 1:4 dilution had better linearity and parallelism closer to 1

Plate format evaluation

Goal: Perform a statistical evaluation of the sample plate




Figure 9. A) Plate layout, B) Design for replicate-to-replicate (Intra-plate), plate-to-plate (Inter-plate), and run-to-run

Summary of Statistical Evaluation for the JM-1 Cell-Based Reporter Assay

Table 1. JM-1 Assay Evaluation

Variance Components Analysis, ANOVA	
Rep-to-Rep (Intra-Plate)	22%
Plate-to-Plate (Inter-Plate)	18%
Run-to-Run	5%
Assay Variability	29%

JM-1 Cell-Based Reporter Assay: Final Plate Format

Goal: Establish the 96-well plate format for the JM-1 NanoLuc assay




Figure 10. Plate format for the JM-1 NanoLuc assay

JM-1 Cell Line: Characteristics of the NanoLuc and plaque potency assays

Table 2. Comparison of the JM-1 and plaque potency assays

Reagent / Item	JM-1	Plaque Potency Assay
Tissue Culture Plate	3 (96-well)	~11 (12-well)
Detection Reagents, volume	Substrate, <2 mL	Coomassie blue, up to 48 mL/plate
Measurement	Automated/ Luminometer	Manual Counting with a light box
Total assay time and data reporting	1 week	1 week

Conclusions

- Developed a Nanoluc reporter cell line, JM-1 for rVSV- Δ G-ZEBOV-GP
- Demonstrated a correlation between the cell-based reporter and plaque potency assays
- Demonstrated JM-1 was stability-indicating for rVSV- Δ G-ZEBOV-GP
- Established an assay for the JM-1 cell line and completed a statistical evaluation

Acknowledgements

Analytical

Melissa Whiteman

Katrina Feller

Kristine Kearns

Amy Bowman

Statistics

Phil Bennett

Formulation

Sara Zimmermann

Jeffrey T. Blue

Robert K. Evans

Akhilesh Bhambhani

Erica Strable

Lynne Isopi

Co-Ops

Vishal V. Ahuja

References

1. Cromeans, TL, Lu, X, Erdman, DD, Humphrey, CD, Hill, VR. 2008. Development of plaque assays for adenoviruses 40 and 41. J. Vir. Methods. 151(1): 140-145.
2. Hall et. al. 2012. Engineered Luciferase Reporter from a Deep Sea Shrimp Utilizing a Novel Imidazopyrazinone Substrate. ACS Chem. Biol. 7: 1848-1857.
3. Patent application US20200181678A1, Cell-Based Reporter Assay for Live Virus Vaccines. B. Meyer, A. Bhambhani, J. Blue, 11-Jun-2020.