




Use of High content imaging to enable Cell-based Potency assay development for ADCs

Madhumati Mukherjee, Ph.D

CASSS Bioassays

Session III-ADCs

April 09, 2025





A biological missile...

1

MOA

2

CMC strategies for potency assays for ADCs

3

Challenges in developing potency assays in the current ADC landscape

4

IncuCyte as a work-horse to enable cell-based potency assay development for ADCs

5


Using the IncuCyte characterization of biological activity

6

Concluding remarks




Primary ADC MoA is payload delivery / warhead-induced cytotoxicity



Standard GMP potency assays for ADCs


mAb intermediate



1. Target-antigen binding

- Cell or non-cell based
- Ensures potency before conjugation
- Common methods: ELISA
- Common readouts: Fluorescence or colorimetric

ADC (DS and DP)




1. Target-antigen binding

- Identical assay as mAb intermediate
- Ensures conjugation does not impact target binding


2. Cytotoxicity assay

- Cell-based
- Common endpoints: reduction in ATP production, membrane integrity
- Common readout: Luminescence or colorimetric or fluorescence




ADCs can possess secondary MoAs


1. Bystander effect




2. Antigen signaling inhibition



3. Fc-effector functions



ADCs can possess secondary MoAs




Control may be required

- Mechanisms should be characterized at early phase.
 - ADCC, CDC, and ADCP
- If activity detected, CMC control is expected (characterization assay).
 - Cell-based functional assays, FcγR binding, and/or N-linked oligosaccharide profile

Frequency and number of commercial approvals is only recently on the rise

13 Commercial Approvals to Date



Next generation ADCs and beyond....

New conjugation chemistries

- Site-specific: Inserted cysteines or unnatural AA's
- Branched linkers: [Anami et al., (2017) Ang Chem Inter 56:733-737]
- Non-covalent conjugations: [Gupta et al., (2019) Nat Biom Engin 3:917-929]



Non-traditional antibody formats

Growing toolbox of warheads

- Greater diversity of MoAs, trigger different pathways to cell death
- Impact of DAR # and DAR distribution to cytotoxicity and effector function
- Wide range of potencies
- Multi-warhead conjugates

More selective target antigens for cancer or tumor microenvironment

Enabling new and accelerated ADC pipelines requires outside-the-box strategies, methodologies, and technologies for CMC bioassay



Using the IncuCyte as
a work-horse to
develop Cytotoxicity
Assays



Live cell imaging: IncuCyte (Sartorius)

Phase and fluorescence high content imager contained within a 37°C incubator

- Continuous monitoring and analysis, allows for early detection and decision making
- Near real-time visualization and automated analysis of living, non-perturbed cells
- Up to 5 fluorescent channels allowing imaging of fluorescently labelled cells (SX5)
- Improves the ability to support newer modalities working under aggressive timelines




Live cell imaging: IncuCyte current uses in our group


In-vitro characterization

- Cell line screening
- Cell seeding density optimization
- Spheroid growth
- Cytotox assay development
- Fc effector function: ADCC, CDC characterization
- mAb internalization studies
- ADC bystander effect- ongoing effort

T-75 Flask





Phase confluency (%) over time




12 hrs

48 hrs


80 hrs



Cytotoxicity assay dev.




CDC assay dev.



mAb internalization

ADC-A RS


isoD102



Finding the correct cell line for the Cytotox assay can be a challenge

Cytotoxicity by Cell-Titer Glo


ADC 1.1 (WARHEAD A, DAR2)



Sensitive

11 (of 11) cell lines


ADC 1.2 (WARHEAD B, DAR8)



- ADCs directed against the same target antigen
- Different warheads with different MOAs and different DAR distribution
- Lack of sensitivity could be due to the warhead's susceptibility to neutral-basic pH
- Screening each cell line using traditional methods in a cytotox assay that takes 4+ days to run would take too long





Example of cell line screening for Cytotoxicity Assays using a 96-well plate in the IncuCyte



Cell culture condition optimization using the IncuCyte

% Confluency vs Time
T-75 cm² flasks
Cell Line A
(no drug)



Days in Culture	Projected Seeding density Flask size (cm ²) (example)	
	T-75 flask	T-175 flask
1	Est. 10.5	Est 24.5
2	6.2	14.5
3	3.1	7.2
4	1.3	3.0

Once a cell line has been selected....

- Label-free confluency analysis
- Perform objective and quantitative analysis of confluency
- Leave cells in the IncuCyte and walk away...
- Extrapolate data to find optimum seeding densities and passaging time
- Gives accurate estimates of doubling rates for cell lines with no prior experience




Cytotoxicity Assay development

Preliminary scouting experiment using the IncuCyte

In this experiment...



- Start to see killing around Day 4
- Plateaus around Day 8
- Can vary depending on warheads, DAR ratios
- Helps to determine optimal assay window
- Live cell analysis vs end-point analysis during assay development work



Spheroid culture increases susceptibility to ADC-induced death

- Growing evidence suggests that culturing cells in 3D better mimics in vivo properties
 - Morphology
 - Genetic epigenetic regulation
 - Metabolism
 - Drug sensitivity and toxicity
- Seeding cells in non-adherent, round-bottom assay plates most GMP friendly way to encourage spheroid formation
- Screening spheroid cultures becomes easy with the IncuCyte
 - Can observe spheroids in 3D
 - Quantify time needed to form spheroids
 - Quantify time needed to kill

ADC- induced cell death





Biological
Characterization
Assessment using
the IncuCyte



Investigating ADC internalization using the IncuCyte-enabling CQA evaluation

Impact of asparagine deamidation to antigen binding and subsequent ADC internalization



ADC-A samples were labeled with 5 μ /g/mL Incucyte® Human Fabfluor-pH Red antibody labeling dye, which is a pH sensitive red-fluorophore that only fluoresces in the low pH environment of a lysosome.

Imaged by Incucyte S3 fluorescent imager for phase contrast and red fluorescence under 20x objective every 20 min.




Fc Effector Function: *In vitro* ADCC assay development using the IncuCyte

Label target cells with Cytolight green marker (live cell cytoplasmic labelling)



Target Cells+ effector cells+ Drug:
Mix labelled target cells with freshly isolated PBMCs (effector cells) at 20:1 E:T ratio and add to assay plate with dose response curve of mAb

IncuCyte: Place 96-well plate into IncuCyte and perform live cell imaging. Real time visualization of target cell killing

ADCC Assay Format



In vitro ADCC activity determination using IncuCyte



Conclusions & Discussion

- **New ADCs hold great promise in the clinic, but pose unique challenges to CMC bioassay**
- **Growing number of warheads require robust toolbox of cytotoxicity methods**
 - Range of potencies and MoAs may require new endpoints and readouts
 - Potency and efficacy of some warheads are not well recapitulated in vivo
 - 3D culturing may overcome in vitro sensitivity challenges
- **Live Cell Imaging using the IncuCyte has greatly improved our efficiency in developing Cell-Based lot release potency assays, characterization assays**
 - Reduces time for screening appropriate cell line that are sensitive to killing by new warheads
 - Employed live cell-imaging to aid in cell line characterization, culture optimization, cytotox assay development work, Fc effector function characterization, CQA assessments



Acknowledgements

AZ-Bioassay Development

Patrick Hussmann

Jing Zhao

Deepa Vedam

Snehal Dighole

Arsala Wallace

Elizabeth Christian

Dylan Weil

Christina Grigoriadou

Scott Umlauf

Sartorius- IncuCyte Team

Lorenzo McMillan

Erin Weddle

CASSS Bioassays Organizers



Thank
you!




BACK UP




Investigating ADC internalization using the IncuCyte-enabling CQA evaluation


A.



B.



C.



ADC-induced death in spheroids


Spheroid-based cytotoxicity assay is...

- Accurate and precise within 50 – 150% RP range
- Stability-indicating
- Demonstrates acceptable performance for qualification and GMP-use


Mechanism of increased sensitivity not well-known, but could be due to...

- Increased antigen surface expression
- Syncing cell replication
- Inclusion of more cell-cell junctions and physiological morphology
- Creation of tumor microenvironment

Linearity



Stability-indicating



Next generation ADCs and beyond....keeps us on our toes!

Evolving conjugation strategies and growing toolbox of linkers and warheads are advancing ADCs through the clinic...


New conjugation chemistries

- Site-specific: Inserted cysteines or unnatural AA's
- Branched linkers: [Anami et al., (2017) Ang Chem Inter 56:733-737]
- Non-covalent conjugations: [Gupta et al., (2019) Nat Biom Engin 3:917-929]

Non-traditional antibody formats

Growing toolbox of warheads

- Greater diversity of MoAs, trigger different pathways to cell death
- Impact of DAR # and DAR distribution to cytotoxicity and effector function
- Wide range of potencies
- Multi-warhead conjugates



Enabling new and accelerated ADC pipelines requires outside-the-box strategies, methodologies, and technologies for CMC bioassay



Confidentiality Notice


This file is private and may contain confidential and proprietary information. If you have received this file in error, please notify us and remove it from your system and note that you must not copy, distribute or take any action in reliance on it. Any unauthorized use or disclosure of the contents of this file is not permitted and may be unlawful. AstraZeneca PLC, 1 Francis Crick Avenue, Cambridge Biomedical Campus, Cambridge, CB2 0AA, UK, T: +44(0)203 749 5000, www.astrazeneca.com




Common methods for ADC lot-release potency testing

Binding

Indirect ELISA/Immunoassay




Competitive ELISA/Immunoassay




Cytotoxicity

CellTiter Glo





Fc Effector Function: In vitro Complement Dependent Cytotoxicity (CDC)



Cells + Serum+ Drug
**Red Fluorescent Marker
for Cell-death (Yo-Yo-3®)**



Using the IncuCyte to optimize CDC activity assay



- Measure fluorescence at multiple time-points within a single plate
- Read plates for days with as little as 5min intervals
- Determine assay window/incubation time

Rituximab Dose-Response at T=7hr Incubation

