



Challenges and Strategies for Potency Assay Development in High-Valency Vaccines



Milica Grozdanovic

Merck & Co., Inc., West Point, PA USA


(28 April 2026)

Overview: Multi-valent and Conjugated Vaccines

- Successful vaccines arise from not only sound scientific strategy, but also from manufacturing and cold-chain supply chain logistics to meet commercial demand
- Optimization of logistics and improved patient adherence is often achieved by the well-established **practice of combining individual vaccines** (e.g., diphtheria/tetanus/pertussis or mumps/measles/rubella)
- Encapsulated bacteria **contain multiple serotypes** – distinct variations within the species (classified based on the structure of their capsular polysaccharide)
- ❖ Conjugated vaccines: **bacterial polysaccharide antigens are linked to carrier proteins** to induce strong immunity in children and adults
- ❖ Revolutionized public health, particularly for *H. influenzae*, *S. pneumoniae*, and *N. meningitidis*





Effect of Pneumococcal vaccines: A game of serotypes



Effect 1: A decrease in the overall rate of pneumococcal pneumonia and associated diseases

Effect 2: Decrease in incidence of vaccine type serotypes but increases in non-vaccine serotypes along with serotype switching.







Regulatory Requirements for Release & Stability

Serotype specific assay



23F



23A



23B



Immunoassays are routinely used for complex multivalent polysaccharide vaccines

- Generation of highly specific mAbs – capture and detection- against structurally very similar polysaccharides (Ps)
- Responsibility to define clinically relevant antibodies, especially when no animal model exists.
- Challenging to obtain since many Ps tend to be highly-cross reactive and can be poorly immunogenic in animals.
- Dedicated teams to handle critical reagents generation, storage, stability

2 X 21 = 42 critical reagents

1

Regulatory Requirements for Release & Stability


Compounding risk for lot release



2

- Measured dose must be 70-130% of the label claim
- 80-120% confidence interval around reportable result
- Challenge:
21 results for release per lot
10 (timepoints) x 21 = 210 results for stability per lot
- 5% failure rate based on assay variability is acceptable for a monovalent product
- If you have a 21 valent product, chances of failing to release a lot raises to **66%**

How we meet the Precision & Accuracy Requirements



Total number of reportable results from ELISA assays in support of Capvaxive development

Capvaxive ELISA by the Numbers

- Number of GMP analysts needed to measure all 21 serotypes in 1 day: **10**
- Time needed to complete assay: **6-7 h**
- Number of reagent solutions prepared: **42**
- Number of microtiter plates: **63**
- Number of analyst interactions with each plate: **12 (756 total)**




Cumulative number of manual analyst steps that need to be executed correctly = 798

Plate randomization & Liquid Handlers




- Plate positioning effects translate into potency bias
- Block randomization mitigates impact of positioning effects on potency

- Liquid handlers for dilution curve plating remove analyst errors and can perform block randomization:
 - Dilution points in the R1-R3 curve zone (A)
 - are prioritized for dispensing in the center of the plate (B)




Block-randomization Results



B

	1	2	3	4	5	6	7	8	9	10	11	12
A	7	7	6	5	5	4	4	5	5	6	7	7
B	6	4	3	3	3	3	3	3	3	3	4	6
C	4	2	2	1	1	1	1	1	1	2	2	4
D	4	2	2	1	1	1	1	1	1	2	2	4
E	4	2	2	1	1	1	1	1	1	2	2	4
F	4	2	2	1	1	1	1	1	1	2	2	4
G	6	4	3	3	3	3	3	3	3	3	4	6
H	7	7	6	5	5	4	4	5	5	6	7	7

“The Potency assay” is an umbrella term for 21 different assays



Examples

Case study #1 Signal to Noise Background

- Blocking buffer containing BSA chosen during development
- Robustness performed for incubation time, BSA concentration, volume/well
- 2 labs reported high Signal/Noise for only one of 21 Serotypes (STs)
- Investigation initially focused on mAbs and samples due to serotype-specific issues
- ❖ New BSA lot caused issues only for one ST!

Case study #2 Low reportable result


- 1 ST reported with results below target
- Investigation initially focused on mAbs and samples due to serotype-specific issues
- mAb lots compared, lab to lab confirmatory testing, plate readers compared
- ❖ Catalog number of solution changed in method accidentally yielded 10X lower BSA conc

(examples are different STs 😊)

Reference Stability: challenges of multi-valent DP


- Concentration of each serotype is often very low
- Very few methods provide type specific information on the multivalent formulation
- Multi-valent formulations can mask drift in individual components
- Variability usually taken from most variable ST

2-Tier Approach - Primary Reference Standard (PRS) and Secondary Reference Standard (SRS)



DS is monovalent and formulated at high concentrations
Monitored using traditional methods (e.g., HPLC, Biophysical methods)
Stored at -70°C and compared against multi-valent Ref Std on stability

Automation: ideal solution for high-valency?



Fully automated system:


- 🦾 Reagents loaded onto the deck of a system equipped with a robotic manipulator arm
- 🦾 384-well multi-channel arm (MCA)
- 🦾 8-channel liquid-handling arm (Tecan EVO).
- 🦾 KiNEDx robotic arm with hoteling stations for storage of plates and tips.
- 🦾 KiNEDx arm transferred plates to the Tecan EVO for each step.
- 🦾 MCA picks up the tips, aspirates and dispenses reagents to plates, washes the tips, and returns them to the box

Challenges:

- 👤 Not QC lab friendly (Vaccine life-cycle duration, testing on import)
- 👤 Automation is more friendly for “parallel” applications and not “serial” (i.e. same assay 20 times, not 20 assays 1 time)
- 👤 Multi-hour setup does not pay off for < 5 samples

Conjugated vaccines: Multi-attribute meets Multi-valent

- ❑ Product is a complex mixture of different molecules, total serotype-specific Ps dose is not enough
- ✓ Need to develop assays that measure only Conjugated-Ps, **Free Ps**, or Free protein



Challenge 1: *Complete removal* of CRM-conjugated Ps prior to measuring Free Ps

- 21 STs have different MW, charges, chemical properties

Challenge 2: Assay has to be very sensitive and serotype specific to measure remaining low Free Ps levels in DP


- < 400 ng/mL for each individual serotype

Magnetic bead immunoprecipitation of CRM-conjugated Ps

- ✓ Our goal: complete depletion of the antigen from the sample.
- ✓ Utilizes 3 clones of anti-CRM mAbs to pull down Conjugated-Ps

Free Ps remains in the supernatant


Free Polysaccharide ELISA



Data from: Grozdanovic, M., Samuel, R., Grau, B. et al. Serotype-specific quantification of residual free polysaccharide in multivalent pneumococcal conjugate vaccines. *Glycoconj J* 41, 47–55 (2024). <https://doi.org/10.1007/s10719-023-10143-6>

Conjugated vaccines: Multi-attribute meets Multi-valent

- ❑ Product is a complex mixture of different molecules, total serotype-specific Ps dose is not enough
- ✓ Need to develop assays that measure only **Conjugated-Ps**, Free Ps, or Free protein



The Multi-attribute ELISA combines 2 separate assays into one assay


Coat: Ps-specific capture mAbs
Load: 11-point DP dilution curve

Detection 1: Anti-CRM mAbs + HRP-IgG


HRP fluorescent substrate

Detection 2: Anti-Ps mAbs + AP-IgG

AP fluorescent substrate





Measures Conjugated-Ps



Measures Total-Ps



Multi-attribute ELISA

A





Left: Representative reference standard curves from serotypes 6A, 11A, 16F, 23B, 31, 35B

B



Right: Average accuracy across 3 dose levels

Multi-attribute ELISA




Order of detection steps matters:

A) The plate was first developed using QB fluorescent substrate followed by 4-MuP

B) The plate was developed first with 4-MuP as substrate, washed, and then developed with QB as substrate.

Conclusions & Acknowledgements



High-valency will compound any issues, is resource intensive, and requires adapting with scale of operations in mind



Would not be possible without support & effort from the whole Merck CMC team

Special thanks to the Capvaxive potency team & friends of the team!!!

Fran Ansbro
Rachelle Samuel
Elizabeth Shenk
Yun Liu
Liming Guan
Gabriella Tamassia

David Holland
Brittney Pedrazzi
Cindy Pauley
Brendan Grau
Thorsten Verch
Jeremy Trausch

Michelle Homsher
Katrina Feller
Jae Roman
Debra Krantz
Kiev Fehnel