

Rapid and High Throughput CE Assays of Enzyme Activity: Liberty Ballroom B&C 13:00-14:30

Lisa Holland,
Lloyd Bwanali, Laura Casto-Boggess, Srikanth Gattu, Jessica
Guerrette, John Lucas
NIH R01GM140560, R01GM114330



Greetings from West Virginia!



Greetings from West Virginia!



Capillary Nanogel Electrophoresis Applied to Glycosylation

Enzyme inhibitor global market \$178B in 2023

Glycosylated substrates are non templated therapeutics requiring sophisticated biotechnology. Complex analytical tools are costly and difficult to use

Self assembled nanogels are multifunctional and improve biomolecule separations at a cost of \$0.76 for the on-board fluid.

Innovation: smart materials enhance functionality at for automated, rapid, and microscale analyses

Findings: Phospholipid nanogels provide **sieving/electrophoresis, enzymatic modification, sequencing, and patterning.**

Goals for Integrated Enzyme Analyses

1. Enzyme Assay: Capillary Electrophoresis is used to quantify enzyme activity and to screen and even quantify the ability of molecules to inhibit the enzyme

Innovation: nanogels pattern nanoliter reactions for automated, rapid, and cost-effective screening

2. Microscale Sequencing: Enzymes are integrated into a capillary to determine glycan (a PTM), sequence and linkage in capillary electrophoresis run

Innovation: enzyme sequencing provides cost-effective identification with no standards or MS

Unique Topics to be Addressed


What makes oligosaccharides difficult to analyze

What is a Nanogel and how does it work?

What makes capillary electrophoresis ideal for integrated enzyme analyses

Why does it matter?

Post-COVID *Significance* of Sialylation (a form of glycosylation) in Scientific Research



Neuraminidase and hemagglutinin function and balance are implicated in virulence. The complexity of these interactions makes it difficult to develop strategies to eradicate influenza A.

Post-COVID *Significance* of Therapeutics and Glycosylation in Scientific Research


Enzyme inhibitors important therapeutics including COVID-19

Glycosylation involved in viral infections
Implicated in binding domain positioning, immune system evasion, cell entry


What makes oligosaccharides difficult to assay or analyze?

What makes oligosaccharides difficult to analyze


monomers




D-Glucose (●)




D-Mannose (●)




D-Galactose (●)



L-Fucose (▼)







N-Acetyl-D-glucosamine (■)



N-Acetyl-D-neuraminic acid (◆)



linkages




The Anatomy of a Glycan

N-glycans (unmodified)


- ◆ sialic acid
- galactose
- *N*-acetylglucosamine
- mannose
- ▶ fucose



Features Monitored



Analysis of Antibody Glycosylation




Unique Topics to be Addressed



**What is a Nanogel and how
does it work?**

**What makes capillary electrophoresis
ideal for reaction-based analyses**


Electrophoretic separations transitioned from slab gel (1959) to capillary gel electrophoresis (1984).



Cross-linked gels in a slab replaced with aqueous solutions or even linear gels in a capillary
+higher voltages
+improved heat dissipation




Electrophoretic separations transitioned from slab gel (1959) to capillary gel electrophoresis (1984).



Cross-linked gels in a slab replaced with aqueous solutions or even linear gels in a capillary
+higher voltages
+improved heat dissipation




Nanogel Electrophoresis (no EOF)



Cunliffe, J.M., N.E. Baryla, and C.A. Lucy, Phospholipid Bilayer Coatings for the Separation of Proteins in Capillary Electrophoresis. *Anal. Chem.*, **2002**. 74: p. 776-783.

White, C.M., R. Luo, S.A. Archer-Hartmann, and L.A. Holland, Electrophoretic Screening of Ligands under Suppressed EOF with an Inert Phospholipid Coating. *Electrophoresis*, 2007. 28(17): p. 3049-55.




Nanogel Electrophoresis (no EOF)



- composed of self assembled lipids
- semi-permanent
- zwitterionic
- can be further modified
- surface is compatible with nanogel
- surface is biocompatible

Multifunctional Nanogels

Fluid Steering, Sequencing and Separations



Phospholipids

5 nm nanodisks

5 nm nanoribbons

Morphology \Rightarrow Viscosity

Thermally Responsive

Temperature Affects viscosity

Liquid at 24 °C

Gel at 26-31°C

Thermo-responsive Nanogels are Easily Introduced and Expelled




Temperature Affects Viscosity

Liquid at 24 °C
Gel at 26-31 °C


Fill/expel “cold”

Lock “warm”


Capillary Electrophoresis Separations of Proteins



Electrophoretic mobility



Sieving




	0% nanogel	vs	25% nanogel
Peptide: REDV (0.57 kDa), pI 4.4,	0.000150 $\text{cm}^2\text{V}^{-1}\text{s}^{-1}$		0.000150 $\text{cm}^2\text{V}^{-1}\text{s}^{-1}$
Alpha-1-antitrypsin (52 kDa), pI 4.6	0.000110 $\text{cm}^2\text{V}^{-1}\text{s}^{-1}$		0.000071 $\text{cm}^2\text{V}^{-1}\text{s}^{-1}$


Crihfield, C. L.; Holland, L. A. **Protein Sieving with Capillary Nanogel Electrophoresis**. submitted.

B.C. Durney, C.L. Crihfield, L.A. Holland, **Capillary electrophoresis applied to DNA: Determining and harnessing sequence and structure to advance bioanalyses (2009–2014)**, *Anal. Bioanal. Chem.*, 407 (2015) 6923-6938.

Neuraminidase Activity/Stability in Lipid



- Pseudo-immobilization
- Non-covalent
- Turnover double
- Activity > 20x
- Lifetime extended months



Gattu, S., C.L. Crihfield, and L.A. Holland, **Microscale Measurements of Michaelis–Menten Constants of Neuraminidase with Nanogel Capillary Electrophoresis for the Determination of the Sialic Acid Linkage.** *Analytical Chemistry*, 2017. 89(1): p. 929-936.

Two Goals for Integrated Enzyme Analyses


Enzyme Assay: Capillary Electrophoresis is used to quantify enzyme activity and to screen and even quantify the ability of molecules to inhibit the enzyme

Innovation: nanogels pattern nanoliter reactions for automated, rapid, and cost-effective screening

Microscale Sequencing: Enzymes are integrated into a capillary to determine glycan (a PTM), sequence and linkage in capillary electrophoresis run

Innovation: enzyme sequencing provides cost-effective identification with no standards or MS

Enzyme Shift Assays are Compatible with Capillary Nanogel Electrophoresis



- Nanoliter volumes of enzyme/lectin (↓cost)
- Pattern in capillary (enables processing)
- Retains enzyme/lectin activity (↓cost/repeatability)
- Replace each run (eliminate carryover)

Unique Topics to be Addressed

What makes capillary electrophoresis ideal for integrated enzymes and reaction-based analyses?


Enzyme Processing with Nanogels

shift in migration time
identifies terminal monomer


no enzyme plug




enzyme plug



Modification with Enzyme Microscale Sequencing




The integration of enzymes in series in a capillary for the hydrolysis and separation of glycans from Trastuzumab.



- Automated process
- nL zones of enzyme
- Pattern in capillary
- Monitoring shifts
- Replace each run


Nanogel Multi-Enzyme Processing to Determine Sialic Acid Linkage



$\% \alpha 2-6$ sialic acid (◇) = $\%$ galactose (●)

- Subtle differences in Sialic acids
- Sialic acids anionic
- Co-cleavage converts 2-3 linkage into an exposed GlcNac
- 2-6 linkage correlates with Gal

Nanogel Multi-Enzyme Processing to Determine Sialic Acid Linkage






- Human AGP (biomarker)
- Sialic Acid Linkage
- 5 zones of enzymes and lectins (nL)
- No Standards

Table 1. Peak area of α 2-6 linked sialic from AGP glycans



Branch	α 2-6 linked sialic acid (◆) area x 10 ³ (%RSD)					Total
	0 ◆	1 ◆	2 ◆	3 ◆	4 ◆	
	25(10)	89(7)	0	—	—	114(12)
	104(5)	336(8)	0	0	—	440(9)
	94(7)	265(3)	0	0	0	359(8)

Nanogel Enzyme Zones Characterize Enzyme Activity




Overview of neuraminidase screening (IC_{50}) and K_I Determinations

In-capillary Enzyme Processing




Improved In-capillary Mixing


Incubate
40s (3 passes)





Incubate
60s (5 passes)



Incubate
100s (7 passes)




Michaelis-Menten Curve for 3'-Sialyllactose



- [E] is constant
- [S] is varied
- Velocity = [P]/(incubation time)

Adapted from Gattu, S., C.L. Crieffield, and L.A. Holland, **Microscale Measurements of Michaelis–Menten Constants of Neuraminidase with Nanogel Capillary Electrophoresis for the Determination of the Sialic Acid Linkage.** *Analytical Chemistry*, 2017. 89(1): p. 929-936.

Enzyme Characterization



2',3',6',8' – Neuraminidase

3'-Sialyllactose

$$K_m = 3.3, V_{max} = 2090 \mu\text{M}/\text{min}$$

6'-Sialyllactose

$$K_m = 2.0, V_{max} = 420 \mu\text{M}/\text{min}$$

3'- Neuraminidase


3'-Sialyllactose

$$K_m = 2.7, V_{max} = 940 \mu\text{M}/\text{min}$$

Adapted from Gattu, S., C.L. Crihfield, and L.A. Holland, **Microscale Measurements of Michaelis–Menten Constants of Neuraminidase with Nanogel Capillary Electrophoresis for the Determination of the Sialic Acid Linkage**. *Analytical Chemistry*, 2017. 89(1): p. 929-936.


Quantification of Inhibitor Activity

$$K_i = \frac{IC_{50}}{\left(\frac{[S]}{K_m}\right) + 1}$$




- K_i sigmoid (dose-response) curve
- Range of Inhibitor Concentrations Used
- Inhibitor Patterned (40 nL volumes)
- Single enzyme preparation (5 nL consumed per run)
- K_i calculated when K_m known of [S] ≪ K_m
- K_i varies with substrate and Inhibitor

Determination of Inhibitor Activity






B. substrate, 3'-sialyllactose conjugated to 2-AB



- Rapid determination of enzymatic rates
- Addition of inhibitor enables study of inhibitor activity against enzyme sialic acid cleavage
- CE fast analysis (5 min), small sample, enzyme (5 nL), and inhibitor volumes

Casto-Boggess, L.D., L.A. Holland, P.A. Lower-Yolar, J.A. Lucas, and J.R. Guerrette, *Microscale Quantification of the Inhibition of Neuraminidase Using Capillary Nanogel Electrophoresis*. *Analytical Chemistry*, 2022.


Universal and Programmable Inhibitor Patterning





- Universally applicable to different molecular inhibitors
- Adding to enzyme preparation consumes excessive enzyme and requires
- Sandwiching enzyme allows robotics to do sample preparation
- Realize full advantage of automation

Casto-Boggess, L.D., L.A. Holland, P.A. Lawer-Yolar, J.A. Lucas, and J.R. Guerrette, *Microscale Quantification of the Inhibition of Neuraminidase Using Capillary Nanogel Electrophoresis*. *Analytical Chemistry*, 2022.

Inhibitor (DANA) and enzyme are mixed in capillary




Inhibitor type 1 (anion) DANA






- 3.3 mM sialyllactose
- DANA-transition state analog
- Pattern capillary 1-cm enzyme zone
- Sandwich between inhibitor zones (7-cm/1-cm)
- No statistical difference if in all 3 zones, versus absent from enzyme

Inhibitor (Siastatin B) and enzyme are mixed in capillary





Inhibitor type 2 zwitterion
Siastatin B



- 3.3 mM sialyllactose
- Siastatin B
- Pattern capillary 1-cm enzyme zone
- Sandwich between inhibitor zones (7-cm/1-cm)
- No statistical difference if in all 3 zones, versus absent from enzyme



Determination of Inhibitor Activity





- Compare conversion with and without bracketing
- No statistical difference between zones 1,2,3 and 1,3
- Equally applicable to neutral and cationic inhibitors

Determination of Inhibitor Activity

DANA






Siastatin B




- Run in series reagent patterning and conversion are programmed
- Each inhibition measurement is normalized against no inhibitor run (% remaining)
- Normalization every third run (10-15 runs total)

Higher Throughput is Possible




Higher Throughput is Possible




- 3 runs, 24 separations of sialyllactose, lactose
- migration time 1%RSD, 2.26 ± 0.03 and 2.38 ± 0.03
- normalized area (55 ± 6) 10% RSD

Determination of Inhibitor Activity

Series vs Parallel





- Run in series reagent patterning and conversion are programmed
- Each inhibition measurement is normalized against no inhibitor run (% remaining)
- Normalization every third run (10-15 runs total)



- Run in parallel reagent patterning and conversion are programmed but are simultaneous
- Each inhibition measurement is normalized against a single no inhibitor run (% remaining)
- Curves are achieved in a single run of 8 parallel electrophoresis separations


Determination of Inhibitor Activity



- Run in parallel reagent patterning and conversion are programmed but are simultaneous

- Each inhibition measurement is normalized against a single no inhibitor run (% remaining)

- Curves are achieved in a single run of 8 parallel electrophoresis separations




Determination of Inhibitor Activity with Fluorescent Substrates




$$K_i = \frac{IC_{50}}{\left(\frac{[S]}{K_m}\right) + 1}$$

APTS labeled 6' sialyllactose (substrate)



APTS labeled lactose (product)






$$K_i = IC_{50}$$

- Change label to APTS
- Reduce [Substrate] to nM regime
- K_m no longer required for K_i determinations

Determination of Inhibitor Activity




DANA Inhibition



- A single run in parallel is performed for a single curve
- K_i values identical to serial analyses


Determination of Inhibitor Activity

Siastatin B Inhibition




- A single run in parallel is performed for a single curve
- K_i values identical to serial analyses

Nanogel Enzyme Processing to Evaluate Enzyme Function



- transferase enzyme
- donor, cofactor, substrate
- activity, preference
- Small volume


Nanogel Enzyme Processing to Evaluate Enzyme Function




- transferase enzyme
- donor, cofactor, substrate
- Polarity mixing
- 30 min incubation

Nanogel Enzyme Processing to Evaluate Enzyme Function


A. Capillary Reaction Zone



B. Electropherograms





C. K_M Curve




- GalT
- K_m 1.26 mM
- Lit report 1.23 mM


Protein Modification



A. Without Enzyme and +/- Lectin



B. With Enzyme and +/- Lectin



Significance: Nanogel electrophoresis increases the functionality and throughput of a microscale separation-based assay

Self assembled phospholipid nanogels provide fluid steering and tunable selectivity of biomolecules at a cost of \$0.76 for the nanogel.

Innovation: Nanoliter reactions are patterned, processed, and quantified repeatedly. Smart materials enhance functionality at a lower cost and higher throughput.

Questions?
Lisa.Holland@mail.wvu.edu

