

Comparability Strategies for Autologous Genetically Modified Hematopoietic Stem Cells (GM-HSC)

Mike Havert

LET'S
RECODE
THE STORY


Disclosures

- Employee and stockholder - bluebird bio

Overview


- Introduce GM-HSC and product pipeline
 - Variable starting materials
 - Analytics under development
- Manufacturing changes and approaches for demonstrating comparability
 - Analytical Equivalence and Split Apheresis
 - Analytics and Manufacturing improvements

bluebird bio products and pipeline




GM-HSC treatment and manufacture overview


HSC collection
*Mobilization
and apheresis*




**Myeloablative
conditioning**



**Drug product
infusion**



**GM-HSCs engraft and
reconstitute blood cells**



Treatment

Manufacture

Collection

Manufacture


Wash and CD34+ Cell
Selection

Culture


Transduction

Cryopreserved Drug
Product


Potential treatments based on genetic modification of stem cells



GM-HSC historical development



GM-HSC manufacture and controls



Allo-HPC/cord blood cell products



CTGTAC Briefing Document

BLA 125397

HemaCord (Hematopoietic Progenitor Cells-Cord)

Table 2: Testing of HPC-C¹

Product Characteristics ²	Testing	Sample (Type and Timing)	Results of Product Testing
Safety	Infectious diseases – Testing Required. (21 CFR 1271.45 through 1271.90)	Maternal peripheral blood obtained within 7 days of cord blood collection – Type and Timing Required. (21 CFR 1271.80(a) and (b))	All tests negative except non-treponemal test for syphilis when confirmatory test is negative. (Cytomegalovirus (CMV) results are recorded.)
	Sterility – Bacterial and fungal cultures – Testing Required. (21 CFR 211.165(b), and 21 CFR 610.12)	HPC-C (pre-cryopreservation) **	CMV – Report No growth
	Hemoglobin	Cord blood*	No homozygous hemoglobinopathy
Purity and Potency ³	Total nucleated cells (TNC)	HPC-C (pre-cryopreservation)	$\geq 5.0 \times 10^8$ TNC ***/ unit HPC-C
	Viable nucleated cells	HPC-C (pre-cryopreservation)	$\geq 85\%$ viable nucleated cells
	Viable CD34+ cells (flow cytometry)	HPC-C (pre-cryopreservation)	$\geq 1.25 \times 10^6$ viable CD34+ cells ****/ unit HPC-C
Identity	Human leukocyte antigen (HLA) Typing	Cord blood	Report
	Confirmatory HLA typing	Attached segment of HPC-C	Confirms initial typing
	Blood Group and Rh Type	Cord blood	Report

¹ Testing, Sample (Type and Timing), and Results are recommended unless specifically noted as required.

² The PHS Act requires a demonstration that the product is safe, pure, and potent.

³ Other purity and potency assays may be considered under the BLA.


* Cord blood = a sample of unmanipulated cord blood. A red cell sample or other cord blood aliquot obtained after volume reduction may be used for testing with appropriately validated reagents or test systems.

** Sample may be obtained before or after addition of the cryoprotectant.

*** Based on 20 kg recipient, a target dose of $\geq 2.5 \times 10^7$ nucleated cells/kg and $\geq 70\%$ post-thaw recovery = 1.7×10^7 nucleated cells/kg.

**** Based on CD34+ cells $\geq 0.25\%$ of TNC prior to freezing.

Variability of CD34+ cell concentration



bluebird bio manufacturing has changed over time

- Facility (CMO)
- Process
- Materials
- Intermediates
- Specification
- Analytical procedures
- Excipients


How do we know the product is still safe?

How do we know the product is still effective?

How do we know the products are comparable?


bluebird bio manufacturing has changed over time

- Facility (CMO)
- Process
- Materials
- Intermediates
- Specification
- Analytical procedures
- Excipients




Technology transfer

- Transfer of manufacturing processes and/or analytical procedures between facilities or laboratories
- Technology transfers take the outputs of process or method development activities and transfer the knowledge to a different location where a process or analytical procedure will be operated.




Technology transfer overview



Regulatory standards for evaluating risk under IND

- CMC changes that impact Safety - 21 CFR 312.42(b)(1)(i)
 - FDA may place a proposed or ongoing Phase 1 investigation on clinical hold if it finds that: Human subjects are or would be exposed to an unreasonable and significant risk of illness or injury
- CMC changes that impact Efficacy - 21 CFR 312.42(b)(4)
 - FDA may place a proposed or ongoing investigation that is not designed to be adequate and well-controlled
 - Applies to all trials, but usually used for phase 3/pivotal



Types of comparability exercises

- Analytical
- Nonclinical / animal model
- Clinical

Analytical comparability considerations


ICH Q5E Comparability of
Biotechnological/Biological Products Subject to
Changes in Their Manufacturing Process

- What samples?
 - Pre-change vs Post-change
 - Full scale vs lab scale manufacturing runs
 - Retrospective or prospective?
 - Lot retains
 - Split apheresis
 - Sample size - how many?
- What are analytics?
 - Lot release (validated?)
 - More than release
 - Characterization tests
 - In process test
 - Stability data
- Statistical assessment regarding equivalence
 - Demonstrating sameness
 - Setting study criteria

Split Apheresis comparability

- Split starting material eliminates person-to-person variability
- Healthy donor
- Test panel (no potency) - can restore gene expression to cells that already have it
 - Total nucleated cells (TNC)
 - %CD34+ cells
 - Viable CD34+ cells


Split Apheresis comparability



* TNC = Total nucleated cells

Analytics and manufacturing improvements


- Develop analytics and bioanalytics
- Optimize manufacturing process in lab
 - Improved product quality
 - Pre-change $\neq$ post-change*
- Technology transfer to CMO
- Update CMC section IND
- Clinical studies



*ICHQ5E - When pre- and post-change products are not comparable, please consult the appropriate regional regulatory authority


Analytics...

CD34+ cells
(drug product)




- Not possible to assess functional activity in final drug product (CD34+ cells) if transgene is not expressed

Analytics...



- Not possible to functional activity in final drug product (CD34+ cells) if transgene is not expressed


Analytics...




- Gene modification corrects arrest at enucleation step in β -thalassemia
- Potency can be measured as a relative increase in % enucleated cells

Image adapted from: hubpages.com/education/General-Considerations-In-Hematology-Blood-Formation-In-Erythropoiesis

Analytics...



Analytics and manufacturing improvements in clinical studies



Summary

- Technology transfer
- Invest in analytics
 - Potency and stability indicating assays
 - Bioanalytics and data science too
 - Less variability is better
- Minimize manufacturing variability (split aph)
- Get clinical experience (introduce changes early)
- Learn from mistakes
 - Small scale models to de-risk

The bluebird flock

Helena Madden
Ilya Shestopalov
Richard Colvin
Leslie Wilder
Matt Hibbert

Cell Analytics

Christine Beaudry
Lauren Beaudin
Brandon Nguyen
Agnes Lin
Steve Duguay

Stem Cell Biology

Gretchen Lewis
Melissa Bonner

Cell Process Development

Erin Meister
Kayla McLeod
Craig Jones
John Pierciey

Leadership & IP

Michael McDonald
Iva Holden
Susan Abu-Absi
Derek Adams