

Regulatory Considerations for Early-Stage Development of Multi-specific Therapeutics

Hiren Parmar, Ph.D.

Biologist/Evaluator, BRDD, Health Canada



Health
Canada

Santé
Canada

Canada

Disclaimer



The views expressed in this presentation are those of the presenter and do not convey official Health Canada policy.

Outline



Review process for clinical trial applications with multi-specifics at Health Canada



Types of multi-specific molecules



Regulatory expectations from a biologics review perspective on what should be included in the submission



Case studies

CMC Review Process for Clinical Trial Applications at BRDD

Office of Regulatory Affairs (ORA)

- Regulatory screening and processing ([Guidance Document For Clinical Trial Sponsors: Clinical Trial Applications - Canada.ca](#))
- Screening information requests for any deficiency





Biotherapeutics Quality Divisions (BQDs)

- Critical assessment of CMC information
 - QOS/IMPD in Module 2
- Quality information requests
- Final decision by BRDD issued within 30 calendar days

Bi-, tri-, and multi-specific molecules

Herrera et. al.,
Trends in Cancer, 2024

Legend



sified / Non classifié

Regulatory Review of Multi-specific Molecules

FDA guidance document

- [Bispecific Antibody Development Programs - Guidance for Industry.pdf](#)

Health Canada uses risk-based approach

- Emphasis on patient safety for early-stage trials
- Harmonize requirements with other regulatory agencies

Review of Early-stage Regular IgGs vs Multi-specifics

- **Cell line development** (construct, host cell)
- **Manufacturing Process** (intermediates, affinity chromatography)
- Raw materials
- **Process-related impurities** (reagents)
- **Product-related impurities** (homodimers, mis-paired species)
- **Characterization** (correct assembly, potency)
- Comparability
- **Specifications** (aggregates, fragments)
- Batch analysis
- Reference standard
- **Stability/in-use stability**

Manufacturing Process

Cell culture

- Intermediates – culture of two mono-specifics
- Single culture with technologies for correct assembly
 - E.g. Knob-in-hole, charge pairing

Purification

- Multi-specific with Fc region
 - Protein A affinity chromatography
- Multi-specific without Fc region
 - Alternative purification strategy
 - E.g. Ion exchange, hydrophobic interaction, Histidine tag-metal affinity

Process-related impurities

Similar to regular IgG

- E.g. host cell DNA, host cell protein, antifoam, insulin, bioburden, endotoxin

Additional impurities

- Specific reagents used during cell culture or purification

Virus clearance

- Not needed if the host is bacterial cells (e.g. BiTEs)

Product-related impurities

Aggregates and fragments

Homodimers

Mis-paired species

Charged species

Characterization

Mostly similar to regular IgG

- Primary, secondary, and higher order structures
- Glycosylation pattern (if applicable)
- Other PTMs as appropriate
 - E.g. oxidized, deamidated
- Size and charge variants

Correct assembly

Homodimers and other mis-paired species

Characterization (cont.)

Biological activity/Potency

- May require multiple assays
- Binding assays may be acceptable for early-stage, but an assay representative of MOA required prior to pivotal clinical trials
 - Capture simultaneous binding
 - Fc functions, if applicable
 - Sensitive to structural changes

Types of assay

- Enzyme linked immunosorbent assay (ELISA)
- Surface plasmon resonance (SPR)
- Reporter gene assay
- Cytotoxicity assay
- Apoptosis assay

Specifications

General tests

- E.g. color, clarity, quantity, osmolality, process-related impurities, microbial safety

Product-related impurity specifications

- Aggregates and fragments
 - Medium molecular weight (MMW) species
- Homodimers and other mis-paired species
- “Report results” may be acceptable at early-stage with justification
 - E.g. Charged species

Potency

Stability

Similar to regular IgG

- Long-term, accelerated, and stressed storage conditions
- Container closure representative of the actual container closure
- Stability data support the proposed shelf life
 - Comparability demonstrated if using supportive batches

Prone to aggregation and/or degradation

- Assays to monitor aggregates, fragments, homodimers, and/or other mis-paired species

In-use stability/compatibility

- Recovery

Submission Examples



Health
Canada

Santé
Canada

Canada

Case study 1: Asymmetric heterodimeric bispecific antibody

- Amino acid modifications to Fc region to promote heterodimer formation
- Gaps identified related to bispecific nature:
 - Homodimer impurities were not monitored
 - Indirectly monitored through SE-UPLC (HMW \leq 5.0%)
 - Removal through CEX during manufacturing
 - Will be monitored as development progresses
 - HC assessment
 - Low levels (HMW 0.3-0.8%)
 - No safety risk from homodimers
 - Design of the bispecific antibody did not allow for formation of a sole T-cell binding molecule
 - Contribute to the mechanism of action

Case study 2: Asymmetric heterodimeric bispecific fusion protein

- Several amino acid substitutions to Fc region including for knob-into-hole
- Gaps identified:
 - Homodimer impurities were not monitored
 - No masses corresponding to homodimers during characterization
 - Hole-into-hole homodimers detected in fermentation supernatant but were removed during purification
 - Knob-into-knob homodimers were not detected
 - SE-UHPLC and non-reduced CE-SDS will be assessed for resolution of mis-paired species as the development progresses
 - Comprehensive risk assessment (biological effects) for each of the homodimers

Case study 2: Asymmetric heterodimeric bispecific fusion protein (cont.)

- Gaps identified:
 - LMW species were controlled as $\leq 15\%$ by non-reduced CE-SDS with no justification
 - Given that this is a fusion protein, potentially include disrupted fragments
 - Characterization of each LMW variant is unavailable at the early-stage (Phase 1 / 2)
 - Strong justification for safety (no potential biological activity)
 - Increase up to $\sim 7\%$ observed during stability of an investigational batch (batch results $< 3\%$)
 - The limit will be re-evaluated as the development progresses

Case study 3: Homodimeric bispecific antibody (IgG1-scFv fusion protein)

- Gaps identified:
 - Surface Plasmon resonance (SPR) based potency assay
 - Measures simultaneous binding, however, not representative of MOA (apoptosis)
 - Advised to develop a functional potency assay prior to pivotal clinical trials

Case Study 4: Homodimeric Bispecific antibody (IgG1-scFv fusion protein)

- Gaps identified:
 - Acceptance criteria for HMW species are 7% at release and 10% on stability
 - Proposed based on the limited data (2 development and 1 clinical batches)
 - Data range 2.6 - 3.2%
 - Drug product stored frozen at -20°C
 - Limits will be re-evaluated once more batch data available

Summary



Multi-specifics pose unique challenges due to their unpredictable nature



Monitoring of product-related impurities at early-stage



Early focus on potency aligned with MOA



Comprehensive clinical trial application package

Pre-submission meetings



For all types of submissions (e.g. pre-CTA, pre-NDS, pre-SNDS)



CMC and clinical teams can be met separately



Can request written feedback only



Follow-up meetings can be organized



brdd.ora@hc-sc.gc.ca

