

Regulatory Considerations for Multispecific Products

Jennifer Swisher, Ph.D.

Supervisory Biologist

Office of Product Quality Assessment III

Office of Pharmaceutical Quality

CDER, FDA

CASSS Strategy Forum, January 26, 2026

Everyone deserves confidence in their *next* dose of medicine. **Pharmaceutical quality** assures the availability, safety, and efficacy of *every* dose.

The views expressed in this presentation do not represent official FDA policy but are the opinions of the presenter.

They are based on experience as a member of OPQ assessment staff and should not be used in place of regulations, published FDA guidance, or discussions with the Agency.

OUTLINE

- Approved Bispecific Antibody-Based Drugs at FDA
- Validation of Multispecific Manufacturing Processes
 - Impurity Control
 - Characterization and Method Validation
- Considerations for Specifications
- Wrap-up

TYPES OF MULTISPECIFICS

- **Bispecific IgGs**
- **Multispecific antibody fragments (joined with linker(s))**
- **IgGs with tandem binding segments added**
- Multispecific fusion proteins
- Multispecific antibody conjugates




FDA APPROVED BISPECIFICS

Blincyto (Blinatumomab)	Amgen	α -CD3, α -CD19	2014
Hemlimbra (Emicizumab)	Genentech/Roche	α -FIXa, α -FX	2017
Rybrevant (Amivantamab)	Janssen Biotech	α -EGFR, α -MET	2021
Kimmtrak (Tebentafusp-tebn)	Immunocore	α -CD3, α -gp100/MHC	2022
Vabysmo (Faricimab-svoa)	Genentech/Roche	α -VEGF-A, α -Ang-2	2022
Tecvayli (Teclistamab-cqyv)	Janssen Biotech	α -CD3, α -BCMA	2022
Lunsumio (Mosunetuzumab-axgb)	Genentech	α -CD3, α -CD20	2022
Epkinly (Epcoritamab-bysp)	Genmab/ Abbvie	α -CD3, α -CD20	2023
Columvi (Glofitamab-gxbm)	Genentech	α -CD3, α -CD20	2023
Talvey (Talquetamab-tgvs)	Janssen Biotech	α -CD3, α -GPC5D	2023
Elrexfio (Elranatamab-bcmm)	Pfizer	α -CD3, α -BCMA	2023
Imdellra (Tarlataamab-dlle)	Amgen	α -CD3, α -DLL3	2024
Ziihera (Zanidatamab-hrii)	Zymeworks/Jazz	Distinct HER2 epitopes	2024
Bizengri (Zenocutuzumab-zbco)	Merus	α -HER2, α -HER3	2024
Lynozyfic (Linvoseltamab-gcpt)	Regeneron	α -CD3, α -BCMA	2025


Bispecific antibody expression

- Can have multiple constructs in one cell line
- Often, engineering of heavy chains favors bispecific pairing



Process Validation

- However, many products utilize multiple cell lines and parallel upstream processes



PROCESS VALIDATION

Multiple upstream processes means:

- Multiple cell banks need to be qualified
- **Multiple seed trains: more raw materials, more time**
- Adequate control needs to be demonstrated for
 - Intermediates (e.g., parent antibodies)
 - Arm exchange (redox) reaction

PROCESS VALIDATION

Regardless of design of upstream, validation of downstream:

- Removal of unpaired and mispaired products- as well as aggregates and fragments- by downstream steps
- Methods should be adequate to demonstrate control both process- and product-related impurities specific to such a process
- Process criticalities based on CQAs; **methods pivotal in control strategy**

CRITICAL QUALITY ATTRIBUTES


Some critical quality attributes to consider for multispecifics:

- **Homodimers** should be characterized to inform criticality (of attribute and of process steps/parameters that affect them)
- Product **aggregates** can be hyperpotent: Synergistic T-cell activation is a potential safety issue
- Products comprised of fragments joined by **gly-ser linkers: O-linked glycans** can be CQAs

Both process and product control may require more precision- and **more assays**- than with a typical mAb

“AHA!”: Another Attribute

- IgG-based sequences used in novel constructs can result in increased accessibility of sequences



- Find reactivity w/ pre-existing “**AHA**”= anti-hinge antibodies
 - Higher levels found in autoimmune populations, cancer patients, e.g.
- Does not necessarily preclude approval, but risk assessment, ADA/NAb assay validation and performance, and discussion with interdisciplinary FDA team are important

CD3 homodimers: Characterization

- Despite accelerated development, solid risk assessment required. Supportive data might include:
 - Nonclinical studies (e.g., cytokine release profile of homodimer)
 - Levels present in clinical batches (and use of these batches)

Additional methods may be required for adequate control of all product-related impurities, given their potential for increased clinical impact for many therapeutics in this class

CD3 homodimers: Clearance/Control

- If assays have either poor or undefined sensitivity for homodimers,
- Cannot assert that no homodimers are present just because levels are \leq LOQ
 - **Levels required for bioactivity can be very low**
 - Does not justify their absence from the control strategy
 - Cannot demonstrate robust clearance/ process control
 - Higher criticality due to uncertainty
 - Engineering reduces but often doesn't fully prevent these and other high-risk impurities

Method performance


➤ Assay characterization

➤ assays must distinguish all product-related impurities

- Each homodimer
- Each unpaired chain
- This may not be the assay for aggregates/smaller fragments

➤ Assay validation

➤ Demonstrate sensitive, reproducible, and accurate measurement of each product-related impurity



Process-Related Impurity Removal

- Processes that require reduction/oxidation reactions or other reaction conditions to ensure proper pairing
- Validation of **process-related impurity clearance** is often straightforward, as small molecules are easily cleared by downstream
 - Must be demonstrated by sensitive, adequately qualified assays

SPECIFICATIONS & Analytical control strategy

- Rigorous control of product-related impurities required vs. typical mAb due to potential for higher risk
 - May require more orthogonal assays due to analytical limitations
 - Narrower acceptance criteria
 - Patient-centric approach for release and stability AC

- Method qualification and precision are paramount due to importance of impurities

Potency Control Strategy

- When simultaneous binding is not obligate
 - Can have separate binding or cell-based assays for separate activities
- When simultaneous binding is required
 - Bioassay(s) should model recruitment of effector function or other unique targeting/ proximity induction
- Surrogate cellular potency assay requires qualification against assay that reflects primary mechanism of action
 - Both may be required in characterization
- **Comparability, pre- and post-approval:**
 - Monitor individual affinities separately: Can be impacted due to cooperative binding from other modules in cell-based assays

Parting thoughts/ Related topics

- Most approved drugs in this class thus far have been granted priority status due to their clinical indications. If so, consider:
 - *Expedited Programs for Serious Conditions/Drugs and Biologics* Guidance for Industry (2014)
 - CDER OPQ MAPP 5015.13:
“Quality Assessment for Products in Expedited Programs”
- Compatibility/in-use stability studies to ensure consistent quality upon preparation and administration
 - Describe companion diluent manufacture and control in BLA
 - Compatibility and in-use stability studies to support labeling

Current Industry Best Practice on in-use Stability and Compatibility Studies for Biological Products



Blümel, M.; Liu, J.; de Jong, I.; Weiser, S.; Fast, J.; Litowski, J.; Shuman, M.; Mehta, S.B.; Amery, L.; Tan, D.C.T. et al. *J. Pharm. Sci.* **2023**, *112*, 2332–2346, <https://doi.org/10.1016/j.xphs.2023.05.002>.




Figure 1: Elements of an in-use study for an IV administration

2024 CASS CMC Strategy Forum is summarized in **Securing Patient Well-Being Best Practices for In-Use Stability and Compatibility Studies**. Liu, Jing, Bluemel, Markus, Swisher, Jennifer, Jules, Felix, de Jong, Isabella, Litowski, Jennifer, Moses, Michael, Weiser, Sarah, Ibrahim, Basma, Clark, Nicholas and Hu, Qingyan (2025) *BioProcess international.*, 23 (4). pp. 19-26.



Acknowledgements

FDA/ CDER/ OPQ

- Marjorie Shapiro (retired)
- Leslie Rivera Rosado
- Kristen Nickens (Amgen)
- Brian Roelofs

