

Combined HILIC and RPLC Peptide Mapping: Two Methods are Better (And Faster!) Than One

Jim Arndt

Senior Scientist, CMC Analytical Development, Characterization



Disclaimer

The opinions expressed in this presentation are those of the presenter and not necessarily those of Teva Pharmaceutical Industries Ltd. or its affiliates (collectively “Teva”). This presentation has been prepared for discussion purposes only. Neither Teva nor any of its employees or representatives make any representation or warranty, express or implied, as to the accuracy or completeness of any information contained herein. The information and examples presented originate from individual experience and may not represent the full scope and/or examples of Teva. Nothing contained within the presentation is, or should be relied upon as, a promise or representation as to the future and Teva expressly disclaims any obligation to update the information if it should change.

Not All Methods are Equal

- **Desire:**

- *More informative data earlier and more often*
- *Inter- and intra-program data trending*


- **Barriers:**

- *Involved and variable sample preparation*
- *Extended chromatography; lengthy data processing*
- *“Force-fitting” the legacy methods*

Solution: *Intelligently design a scalable, selective method for complete characterization*

Unifying and Platforming

Peptide mapping using LC-MS/MS is the **only assay** that provides information on a **per-residue basis**




Embrace Change: Build a superior platform method from the ground up.

Three key aspects:

1. *Sample Preparation*
2. *Sample Analysis Method*
3. *Data Processing and Reporting*


Acid-labile Surfactant for Protein Denaturation



- Surfactant reconstituted in 50 mM ammonium bicarbonate (AmBic)
- All reduction/alkylation/digestion is performed in one pot
 - No risk of sample loss
- After acidification, samples are centrifuged and decanted to remove water-insoluble product
- **Removes the need to perform buffer exchange step**

Image: Waters Application Note 720003102, June 2009

Full Teva-mAb Sequence Covered in Surfactant Method




Legacy Preparation

Surfactant

Proteins	Number of MS Peaks	MS Peak Area	Sequence Coverage	Abundance (mol)
1:LC	117	60.6%	100.0%	60.41%
2:HC	181	39.4%	94.2%	39.59%
Unidentified	0	0.0%		

Small, hydrophilic peptides not retained

Scalability



- Surfactant platform method was used to characterize **low concentration** drug product
 - 136 pmol (20 µg) starting material used
- Obtained full sequence coverage and maintained ability to quantify critical PTMs




Next Steps

- A simplified, platform, scalable digestion protocol has been developed
 - Only 1/3 of the story. . .



- High-resolution mass spectrometry provides additional mode of “separation” with near infinite signal-to-noise
- *Opportunity to significantly shorten analytical method without sacrificing data quality*

Leveraging HRMS to Increase Throughput



Achieving Complete Characterization

- Tryptic peptide maps of IgGs contain both very hydrophobic and very hydrophilic sequences
- 100% sequence coverage is desirable, but not always feasible by a single technique


Proteins	Number of MS Peaks	MS Peak Area	Sequence Coverage	Abundance (mol)
1:LC	117	60.6%	100.0%	60.41%
2:HC	181	39.4%	94.2%	39.59%
Unidentified	0	0.0%		

- Some CQAs, such as glycoforms and deamidations are difficult to resolve, so others have chosen longer reversed-phase gradients or alternative enzymatic digestion

What if we could do two gradients that employ different phase characteristics to provide complete characterization in a shorter amount of time?


Reversed-Phase and HILIC as Complementary Techniques

HILIC



- Polar (amide) stationary phase
- Separation achieved by hydrophilic partitioning into surface water layer
- Good for retention of *polar* analytes
- Mobile Phases:
 - A: water, 10 mM ammonium formate, 0.05% formic acid, pH 4
 - B: acetonitrile, 0.1% formic acid



Reversed Phase




- Non-polar (C18) stationary phase
- Separation achieved by partition into stationary phase with organic elution
- Good for retention of *non-polar* analytes
- Mobile Phases:
 - A: water, 0.1% formic acid
 - B: acetonitrile, 0.1% formic acid

Image: Restek, Inc. How to Avoid Common Problems with HILIC Methods. Lit. Cat.# GNAR2716-UNV. 2017

Science Plumbing!



Two Rapid, Complementary Gradients Provide Full Coverage in Half the Time



Teva mAb
Legacy Method
10 µg injection
128 min. run time
95.8% Seq. Covg.

Teva mAb
New RPLC
5 µg injection
32 min. run time
100% Coverage


Teva mAb
New HILIC
5 µg injection
37 min. run time
98% Coverage

**Combined method covers
100% of the Teva mAb
sequence in ~54% of the time
as the previous method**


Ensuring Continuity Between Methods

- Method bridging study
 - N = 6 preparations of Teva-mAb reference standard prepared by the legacy method and surfactant method with combined RPLC/HILIC analysis
 - Samples prepared using the same reagents, where possible, and analyzed on the same instrument within 24 hours of each other
 - All regularly-reported PTMs are covered by combined RPLC/HILIC approach with similar results to legacy average

Teva-mAb Fc Met Oxidation




Teva-mAb Fc Asn Deamidation



**Error bars represent ± 3 s.d. around mean*

New Method Maintains Ability to Monitor Isomerization CQA in CDR


CDR **Asp**-Gly Peptide; 22 Residues



D



Iso-D




Deamidation Resolution in Fast HILIC

RPLC-MS/MS



Fc **Asn**-Gly Peptide; 16 Residues



HILIC-MS/MS

★ Possible D isomerization


Asp/Iso-Asp-Gly

HILIC MS/MS provides more informative data in < 7 minutes

Asn(-NH₃)-Gly

Heat Stress – Fc Deamidation Uncovered in Small Peptide by HILIC


Asn-Lys Peptide; 6 residues




Teva-mAb; 2 Wk, 37 °C

Teva-mAb; 8 Wk, 37 °C

Asp-Lys, Modified Peptide




Asn-Lys, Native Peptide



Asn-Lys, Native Peptide

Asp-Lys, Modified Peptide

Fc Deamidation Peptide, Heat Stress



Glycopeptide Resolution in HILIC-MS/MS



HILIC-MS/MS can be used to routinely monitor glycosylation profile

A1G0F

A2G0F

A2G1F

M5

M6

M7

M8

M9

High-mannose species detected and resolved



Shortened Gradient with Targeted Analysis


- Late-stage products have well-defined list of PTMs
 - Control chart

	HC PyroE	LC PyroE	HC Fc Deamidation (%)	HC Fc Aglycosylation (%)	HC CDR Ox	HC Fc Ox
Average	1.06	0.27	4.52	0.61	4.86	1.63
SD	0.09	0.03	0.41	0.07	2.04	0.39
%RSD	8.1%	12.8%	9.1%	11.4%	42.1%	24.1%
3SD	0.3	0.1	1.2	0.2	6.1	1.2
Upper Limit (Avg+3SD)	1.3	0.4	5.8	0.8	11.0	2.8
Lower Limit (Avg-3SD)	0.8	0.2	3.3	0.4	-1.3	0.4


- Use of mAbs with conserved sequences lead to the same peptides
- *Build a fully targeted data processing method that can be used across all programs*
 - *Screening assay for early-phase programs*
 - *Routine monitoring assay for late-phase programs*

Inter- and Intra-Program Trending

Fc Met Ox %Mod by Test no., Inter-program



Average Fc Met Ox Value, Inter-Program



Where do we go from here?

- More informative data, *globally*
- Peptide mapping is considered a “Multi-Attribute Method” (MAM)
 - Has the potential to replace multiple assays in a single experiment
 - Most specific assay
- Barrier to adoption has been sample preparation, data analysis, and instrument maintenance
 - One-pot surfactant method significantly decreases sample preparation difficulty
 - Fast gradients and targeted data processing methods enable near hands-off data processing
 - Method is completely scalable to meet varying sample concentration needs
 - Still room for optimization

This work has laid the foundation to take characterization assays outside of the AD lab and into the hands of our global colleagues on a single, harmonized platform

Acknowledgements and Thanks

- AD Characterization Team
 - Mark Platt
 - Kristin Phipps
 - Ozlem Onder
- AD Forced Degradation Group
 - Ervinas Gaidamauskas
 - Rachel Kimble



Thank you.

