

Elucidation, Characterization and Monitoring of a Unique Tyrosine Sulfation Post-translational Modification During Bispecific Antibody Process Development and Scaleup


David Mahon, Nandakumar Madayiputhiya, Ming Gu, Guodong Chen, Li Tao
Analytical Development and Attribute Sciences
Biologics Development

Sept 6th, 2023


CASSS-Mass Spectrometry

Molecular Design, Developability and Biotransformation

Protein Characterization Plays an Integral Role in the Development of Biologics Products



MS Role in Protein Characterization



Presentation Outline

✓ Tyrosine Sulfation Characterization

- Tyrosine Sulfation Introduction and Molecule Background Information
- Identification Site of Sulfation
 - Multi-Enzymatic Proteolytic Digestion Approach
 - Synthetic Peptides
 - EAD Fragmentation


✓ LC-MS for Cell Line Development & Process Development Support

- High-throughput Sulfation LC-MS Intact Mass Method Development
- LC-MS Support of CLD and Process Development Efforts to Reduce Sulfation During Protein Therapeutic Process Development

Sulfo-Tyrosine- Background

Recent Literature

- Tyrosines are sulfated by tyrosylprotein sulfotransferase (TPST) in the Golgi
- TPST motifs are associated with adjacent E/D residues
- AEX purification step removed sulfo-tyrosine species (lower yield).
- Sodium chlorate at 16 mM inhibited tyrosine sulfation by more than 50%- no impact on antibody titer or quality.



Seibert and Sakmar. Pep.Sci. 2008, 90, 459-477

Zhao et al. Mabs, 2017, 9, 985-995


Tyshchuk et al. Mabs, 2019, 11, 1219-1232

Liu et al. Biotech J. 2021



Discovery Data

Phosphatase & Sulfatase Treatment Confirms Sulfo-Tyrosine

2+1 Bispecific




Intact Mass Analysis




Identifying Site of Sulfation- Challenges for Quantitation and Localization Using MS

- Sulfotyrosine is highly labile: MS Source Parameters Can Induce Loss of SO₃




- Complete Loss of sulfate SO₃ (80 Da) Prior to Fragmentation-CID or HCD.
- Impossible to Identify Site if ≥ 2 Y Present on the Peptide



Summary of Multi-Enzymatic Proteolytic Digestion Approach to Identify Site of Tyrosine Sulfation

2+1 Bispecific Antibody



(1) Trypsin digestion identified a +80Da modification on tryptic peptide:

R.XXXXDY¹Y²Y³DXXXXXXXXY⁴Y⁵Y⁶XXX
XXXXXXXXXXXXXXXXXXK.X



Identification with ETD

(1)

(2)

(3)

(4)

(2) Asp-N digestion identified a +80Da modification on HC CDR peptide:

DY¹Y²Y³.D


(3) Chymotrypsin digestion identified a +80Da modification on HC CDR peptides:

Y.Y³DXXXX.X and Y.Y²Y³DXXXX.X

(4) Synthetic peptides confirmed Y³ on HC CDR region as site of sulfation

Synthetic Peptides Confirmed Tyrosine Sulfation on Y³


- XIC of Chymotrypsin Peptides Y³-X and Y²-X +80Da





- Sulfated Peptide in Sample Digest Elute at the Same Retention Time as Y³ Sulfated Synthetic Peptides

Orthogonal Electron Activated Dissociation (EAD) Confirmation: Site of Sulfation on Y³

- Alternative MS/MS fragmentation such as Electron-transfer dissociation (ETD) and Electron-capture dissociation (ECD) are milder fragmentation techniques compared to CID and HCD.
- ECD/EAD generated signature ion (C₈-C₁₀) to locate sulfation on Y³ of tryptic peptide XXXXDYY³DXXXXXXXXXXYYXXXXXXXXXXXXXXXXXXXXXK.X.



Tyrosine Sulfation Impacts Binding of Bispecific to Target




	1% sY change corresponds to change in potency of y% (slope):
sY1	-0.312
sY2	-0.752

- Tyrosine Sulfation Determined to be a CQA


LC-MS Cell Line Development & Process Development Support

High-throughput Sulfation LC-MS Intact Mass Method Development for CLD and PD Support

A: Glycosylated Intact





B: De-glycosylated Intact



Sample Preparation: PNGaseF reduces sample complexity and possible interference from glycans


High-throughput Sulfation LC-MS Intact Mass Method Development: Comparison of Intact MS Approaches

Sample Preparation	Type of Analysis	Reason
 De-Glycosylated	Reduced	<ul style="list-style-type: none"> • Better resolution. • No interference from glycosylation • <u>More sensitive to MS conditions.</u> • Additional data analysis Vs intact-MS
 De-Glycosylated (PNGaseF, 2hr, optimized)	Non-reduced	<ul style="list-style-type: none"> • Better resolution • No interference for sY1, Interference from glycation for sY2 quantification. • Less sensitive to MS conditions

MS Parameters Impact on sY1 quantification		
MS Parameters	De-glycosylated Non reduced	De-glycosylated Reduced, HC1 & HC2
Cone Voltage	Low	High
Capillary Voltage	Low	Low
Capillary Temperature	Low	Low





Impact of Cone Voltage on sY1



Relative Quantitation of Sulfo-Tyrosine During Clone Selection



Goal: Selected a clone with relative percentage lower than reference lot.



- LC-MS based intact mass method to quantify Tyrosine sulfation trends during CLD



Previous Study: Chemical Inhibitors to Reduce Tyrosine Sulfation

Conclusion: Sodium chlorate at 16 mM inhibited tyrosine sulfation >50% with no major impact on antibody titer or quality





Liu, R., Zhang, Y., Kumar, A., Huhn, S., Hullinger, L., & Du, Z. (2021).
Modulating tyrosine sulfation of recombinant antibodies in CHO cell culture by host selection and sodium chlorate supplementation. *Biotechnol. J*

Sodium Chlorate Sulfation Inhibition Study



- Minimal impact on the reduction of sulfation with the addition of 30mM sodium chlorate.
- Sodium chlorate impacted titer.

CLD Re-engineered Cell Line Prevented Tyrosine Sulfation



Summary

- ✓ Tyrosine sulfation presented multiple analytical challenges for identification and quantitation.
 - ✓ Identification Site of Sulfation:
 - Multi-Enzymatic Proteolytic Digestion Approach
 - Synthetic Peptides
 - EAD Fragmentation
- ✓ Structure/Function Characterization of Sulfated Species Including Assessment of their Impact on Potency
- ✓ LC-MS Intact Mass Analysis to Support CLD and PD Efforts to Reduce Sulfation During Protein Therapeutic Process Development.
 - Sodium Chlorate and CLD Study
- ✓ Establishment of Fundamental Guidelines for Identification and Quantification of Tyrosine Sulfation PTM in Therapeutic Proteins

Acknowledgements

BPC CoE

Nandakumar Madayiputhiya

Ming Gu

Guodong Chen

Naresh Chennamsetty

Ran Furman

Thomas Slaney

Li Tao

Bioassay CoE

Patrick Wong

Zhijie Cheng

Darren Kamikura

Tara Stauffer

Biologics Analytical Development

Hao Luo

Ming Zeng

CLD

Jacquelyn Schulman

Dan Huang

Global Upstream and CLD

Evan Schlaich

Rachel Egan

Lakshmi Kandari

DPD

Haiying Bao

Ji Zheng

GPAS

Jessica Gruskos

Analytical Organization Development

Joseph Huang

Hui Wei

DPAS

Christopher Lietz

Sarah Franc

Younhee Cho

Steven Wang

Sally Karbo

Jon Cordia

Doug Banks

PAS

Ekaterina Deyanova

SCIEX

Zhengwei Chen

