

From Blur to Clarity: Definition of Particle Visibility Threshold in Parenteral Drug Products

Felix Nikels, Boehringer Ingelheim

Atanas Koulov, Clear Solutions Laboratories

On behalf of the EFPIA Global Visible Particles Workstream



Introduction



Control of visible particles is an important aspect of the control strategy for injectable products to reduce potential risks to patients



A major technical hurdle for visual inspection stems from different degrees of visibility of particles, rendering particle detection a probabilistic process



To gather new data on particle visibility and shed light on this decade-old challenge, a multi-company blinded visual inspection study was conducted by 10 independent global companies



A major goal of the study was visibility assessment of several particle types of different sizes in small volume vials, as a challenging configuration for visual inspection.


Study Design

Protein surrogate (SU-8), fibers, glass, metal, rubber, PS beads, controls
1 or 5 per container, 50 – 400 μm , 1 ml in 2ml glass vials


Randomized and blinded particle test kit
shipped sequentially to individual companies

Kit was individually inspected according to
EP 02.9.20 by at least 3 qualified analysts per company

Readout: Particle presence (y/n) + additional observations
→ 5460 individual results from 39 analysts




Probability Of Detection (PoD) As Function Of Particle Type And Size



- General trend: presence of multiple particles resulted in higher PoD than single particle per unit
- Detectability of particles smaller than 250 µm has been shown to be unreliable in this study
- Certain particle types, (e.g., PS beads) or multiple particles per unit, have a higher PoD
- PoD >70% reached/exceeded for some particle types between 200 µm and 250 µm, but PoD of fibers did not reach 70% even at 400 µm

Heatmap Overview Of The Multi-Company Threshold Study Outcome



- Average PoD values across companies ranged 5 to 92%, only certain particle types detected with PoD >70%
- PoD increases with increasing number and size, but only for certain particle types
- Counting accuracy was found to be between 0 to 79% (with only 2 instances above 70%)
- Classification accuracy was overall low (with only PS beads and rubber reaching medium values)

Summary



Setting a generic low threshold for visibility of particles (50 μm , 100 μm or any other exact size) is not practical, and likely not achievable, as it will strongly depend on a variety of factors



Under the conditions of the study classification of particle identity and counting of particles by visual inspection alone is very challenging



Study findings underscore the necessity of robust procedures for analyst training and qualification, setting realistic expectations on size-based visibility limits, and harmonization of guidelines globally



Multiple particle control strategy layers are required to provide assurance for a consistent quality of manufactured parenteral products

All graphs were derived from Mazaheri et al, 2023, J Pharm Sci (DOI: [10.1016/j.xphs.2023.10.002](https://doi.org/10.1016/j.xphs.2023.10.002))

Meet The Team

EFPIA Global Visible Particles Workstream



Christof Finkler



Boehringer
Ingelheim

Felix Nikels

Genentech

Miguel Saggiu



Klaus Wuchner

AstraZeneca

Kristen Gonzalez



Bristol Myers Squibb

Tapan Das

Lonza

Pascal Chalus

AMGEN

Jonny Parsons



Maryam Mazaheri

clearsolutions
LABORATORIES

Atanas Koulov



Global Biotech Experts, LLC
Experienced - Engaged - Efficient - Effective

Patricia Cash



Sophia Levitskaya-
Seaman



John Case

Mazaheri et al, 2023, J Pharm Sci (DOI: [10.1016/j.xphs.2023.10.002](https://doi.org/10.1016/j.xphs.2023.10.002))

Special Thanks To All Analysts & Lab Supervisors For Their Essential Contributions!

- **MacroGenics:** Pablo Dominguez, Amaryllis Kanyabwero, Nicole Utano
- **FibroGen:** Julia McNeal, Amy Tse, Ron Merida
- **Bristol Myers Squibb:** Rekha Vasudev, Keyur Sony, Sam Mathew, Eric Payne
- **AstraZeneca:** Sonia Dragulin-Otto, Robert Custer, Arian Torabi, Brianna Stuckey, Marena Martinez-Correa, Katherine Elliott, Colton Philips, Sam Pierre
- **Amgen:** Sean O'Rourke, Luke Finnegan, Johnson Ewulonu, Barry Donohoe
- **Boehringer Ingelheim:** Thomas Ehrmann, Silke Hoevel, Cornelia Gapp, Narine Gaspajan, Olga Carman, Kristina Koepp, Stefan Strecker
- **Janssen:** Debora Bonvin, Luca Schmidle, Auel Tschamer, Karen Neuhold, Susan Claussner, Anna Infortuna, Bozena Pech, Shermine Ljatifi
- **Genentech:** Maggie Yu, Leiza Lopez, Judith Garcia
- **Roche:** Margarete Jehmlich, Juliette Nollet, Ilka Hetzger, Silvia Malecha, Miriam Schott, Saras Dondu, Franka Hosan, Thomas Raetz, Silvia Peter, Larissa Schwab
- **Lonza:** Hatim Chenkaoui, Laura Garcia Rodriguez, Filip Maciej Fedorowicz
- **NIST:** Michael Carrier



From Blur to Clarity:

Definition of Particle Visibility Threshold in Parenteral Drug Products

Atanas Koulov, PhD

Chief Scientific Officer

atanas.koulov@clearsolutions-labs.com




Can we define a universal particle size “visibility” threshold?

“Visibility” limit from a physiological perspective

Physiology of human visual perception: detection of objects in the retina of the human eye


Is based on **2D projections of 3D objects on the retina**

Goldstein, B. E.
Sensation and Perception, 2009




Hypothesis:

2D, area-based particle measurements (e.g. ECD) may be better suited to define particle “visibility” than uni-dimensional measurements




- Inspection panel:
 - 10 qualified operators, 3 inspection rounds each
- Defect test set:
 - 2 different containers - 2R and 20R vial
 - Defect and defect-free units
 - 5 different particle types (metal, glass, rubber, white fibers, black fibers)
 - Gradient of particle sizes in each particle class
- Inspection conditions - according to USP <790>/ Ph. Eur 2.9.20:
 - B&W inspection board with a non-glare white panel base
 - illumination 2'000 - 3'750 lux;
 - 5s/ 5s inspection duration

1D size parameter (Ferret max)




2D size parameter (ECD)



Fedorowicz et al, (in review)

- If 1D parameter (Ferret max) is used - the “visibility” limit varies per particle type
- If 2D/ area-based parameter (ECD) is used - **PoD of 70% (visibility limit) is ca. 200 μm for all particles**

Results (2R vial 20R vial)



Fedorowicz et al, (in review)

Container dimensions affect the particle “visibility limit”:

Larger container sizes results in a larger particle threshold of detectability (“visibility”)



- Using the area-based size-parameter ECD allows for definition of a single “visibility” limit for particle defects of various types and morphologies. This finding presents exciting opportunities for standardization.
- The “visibility” threshold is not absolute. Changes in container dimensions and other product attributes (data not shown) result in shifts of the “visibility” limit.



 expertise. insight. agility.

Thank you!



Filip Fedorowicz

Andreas Zerr

Roman Mathaes

Matthias Eisele

Swen Maas